

FORTH

USER GUIDE

to programs for interactive image processing
and microdensitometry on the RGO PDS.

JONATHAN MCDOWELL

KEITH TAYLOR

RGO 26/07/78

INTRODUCTION

[0.1] Conventions:

- i) User typed FORTH instructions in capitals.
- ii) Terminal response in italics
- iii) Comments and descriptions in lower case.
- iv) In user typed FORTH instructions, lower case symbols such as: n, m, x, y stand for user typed numbers and are defined in accompanying comments. Numbers are to be assumed to be in 16-bit format (integers, unpunctuated) unless stated otherwise.
- v) PDS coordinate directions defined by x and y.
- vi) Size of 2-D data arrays on disk or tape defined by ρ (for Rows) and κ (for columns).
- vii) Pixel identification within those arrays given by i, specifying row number and j, specifying pixel in row (column number).
- viii) Some parameters are specified by 2 components. In this case they will be identified using vector notation in the comments. For example, WN is a parameter representing the dimensions of a 2-D array or scan. Its components are 2, 16-bit numbers referenced by WN (κ - direction) WN $l+$ (ρ - direction).

[0.2] Instructions:

FORTH instructions are single words, delimited by spaces. The words will be executed on typing a carriage return. (represented here by (cr)). On completion of execution, if execution has not caused any errors, output will be *OK*. Several instructions may be typed on line, separated by spaces. They will be executed in turn. FORTH uses a push-down stack for number operations. If an instruction is encountered which is not recognised, the word is echoed back followed by a question mark, and execution is interrupted.

[0.3] Number entry:

FORTH uses single-precision (16-bit) and double precision (32-bit) signed integers. 16-bit integers are in the range $-32768 \leq n \leq 32767$. These are used for most purposes but sometimes, 32-bit numbers are needed, (for example, PDS co-ordinates). 16-bit integers are entered with no punctuation except a leading minus sign where applicable; ie the same as FORTRAN integers (but there must be no leading plus sign). 32-bit integers are entered by including punctuation either , . - or / . Any number which includes punctuation (except for a single leading minus sign which means a negative number for both 16- and 32-bit numbers) will be treated as a 32-bit number, but it will be treated as an integer, ie position of the punctuation will be ignored, hence the numbers 12.34 1.234 12,34 are all taken to be the same number, a double-precision 1234. To avoid confusion, it is probably simplest just to terminate with a point (eg 1234.)

Jonathan C. McDowell

[0.4] Basic operations:

Arithmetic is reverse polish; integer only.
Typing a number puts it on the stack.
Typing . prints out the top of the stack.
Typing (number) (address) ! stores (number) in (address).
Constants To put value of constant, R on stack just type;
R

n ' R! stores n in R
Variables Typing the name of a variable, DX puts address of DX on stack.

To put value of DX on stack type;
DX @

n DX ! stores n in DX
The command ? is equivalent to @ . hence it prints out value in address currently on stack; eg, DX ? prints current value of DX.
To set number base, type DECIMAL, OCTAL or HEX.

For more detailed explanation of subjects in [0.2] to [0.4], see FORTH course notes , FORTH Reference Manual.

[0.5] Note on FORTH words:

The FORTH compiler only stores the first three letters of each word and its length thus, FILE is treated the same way as FILL; FORTH cannot tell the difference. FILE1, FILE2, FILE3, FILE4 are treated the same, but FILEAB is not (different length). This is useful when a long word is spelt incorrectly but can lead to confusion.
Also note: CENTRE = CENTER, DISK = DISC etc. This applies to 'words', ie FORTH commands, but not to alphanumeric strings. These strings are usually used as identifiers or file names. They are only valid if used after a command requiring such a string; such commands include FIND, WRITE, INTER, ENTER. In this documentation, notation of the form (name) or (filename) will refer to one of these strings (eg; FIND (filename) might be FIND GALAXY). These strings should not have any spaces in them, except where mentioned in the documentation.

[0.6] Interrupt:

"Control C" acts as an interrupt. This may not work when tape is being used.

[0.7] PDS control:

For user to move PDS stage, the SELECT (blue and white) button must be on 'manual', for the computer to move it, it must be on 'auto' (eg when doing a scan). Remember to turn control over when using a command which requires stage movement, as the documentation does not mention this. Control can be turned over before or after the relevant command has been sent.

[0.8] Images:

The same basic instruction, PASS, is used to transfer scan data from PDS to tape, or from disk to terminal display, etc, ie from one representation to another. The format of the instruction is (unit1) SOURCE (unit2) DESTINATION PASS to transfer data from (unit1) to (unit2); the format (location, array size, scan sample interval, etc.) is determined by a set of parameters; each unit has its own set. The whole instruction may not always be typed explicitly, eg in plotting routines the word PLOT may be defined as SOURCE SCREEN DESTINATION PASS, and so DISK PLOT will transfer the data from the disk to the screen in graphic form. Each unit name is called an 'image'. A list of 'images' follows.

Images

PLATE The PDS stage; PLATE SOURCE (unit2) DESTINATION PASS does a scan
DISK A file on disk. The word DISK is used implicitly in the instruction
 FIND (filename) which also selects a specific file.
TAPE A file on tape
TERMINAL The terminal; for digital output
SCREEN The Tektronix terminal; for graphic output

Also:
 INPUT)
 OUTPUT)
 OPERAND)

Low-level functional images used by PASS

Parameters
LOCATION

(these can be treated as 'variables')
 Disk block number for DISK image; not used for other images.
 Value is $500 \leq n \leq 2435$.
 Size of scan in pixels; CN contains number of pixels in a row
 CN l+ contains number of rows in the scan
 Sample interval in microns:
 DX contains sample interval in X
 DX l+ contains sample interval in Y
 Origin of scan in PDS co-ordinates. This is the first pixel in the array, designated pixel 0, row 0, (for an nxn array pixels and rows are numbered 0 through n-1.) C0 contains 2 double-precision numbers: C0 contains the X co-ordinate, C0 2 + contains the Y co-ordinate. A 'window' or subset of the array (on disk, tape, etc) may be defined. The origin (first pixel of first row) in pixels relative to the origin of complete array is contained in W0. W0 contains origin in pixels, W0 l+ contains origin row number. ie "i", "j" of origin respectively, see [0.1]. W0 and W0 l+ are always positive.

CN

DX

C0

W0

WN Size of window in pixels; WN contains number of pixels/row of window,
WN l+ contains number of rows in window.
KIND Specified direction of scan, ie scan in X or scan in y.

Current image

The parameters listed above are not ordinary variables because each of the 8 images mentioned above has its own set of parameters. For example, typing DX will reference the DX of only one of these images, called the "current image". An image is made the current image by typing its name (or using a word which itself used an image name). For instance, TAPE DX ? will print out the current DX value for TAPE, while PLATE DX ? will print out that for PLATE, which may be different. Words which change the current image are:

FIND makes DISK the current image
TWIDE makes TAPE the current image
PASS makes OUTPUT the current image
PLOT makes OUTPUT the current image after setting OUTPUT to be equal to SCREEN
SOURCE Copies parameters from current image to INPUT and leaves the latter as current image.
DESTINATION copies parameters from current image to OUTPUT and leaves the latter as current image.
SECOND copies parameters from current image to OPERAND and leaves the latter as current image.

USING FORTH

[1] Bootstrap, Loading

[1.1] Follow 'starting up' instructions in black PDS manual up to point of pressing HALT switch on console; use FORTH SYSTEM DISK 1

[1.2] Disk Bootstrap

Set up correct address on data switches, as in PDS manual, and as indicated on console.

Press console switches: HALT (ie depress HALT/ENABLE)

LOAD ADRS

ENABLE (ie raise HALT/ENABLE)

START

If
(Otherwise

HELLO ? Go to [1.4]
continue on to [1.3])

[1.3] Bootstrap Fails

Hit "Return" key

OK Go to [1.3.1]

Try reloading the disk; see 'closing down' and then 'starting up' again in black PDS manual,

If
Otherwise

[1.3.1] Block Buffers Corrupted

33350 4010 ERASE (should come back OK)

Return to [1.2]

[1.4] Loading

Response is the system date and time, these may be corrected if desired.

JCM LOAD 20/07/78 14:12 OK

Response is the system date and time, these may be corrected if desired.

Sets date to be 1978 Sep 5

Sets time to be 18:07; GMT and BST

are alternatives to UT and are

identical, computer does not know

about time zones.

If you do not wish to correct date and time (this is completely optional) go to [1.5]
Otherwise, eg

5/09/78 NOW OK

18:07 UT OK

To print out system data, time TODAY @ .DATE 5/09/78 OK
 e.TIME 28:08 OK

[1.5] If required, load magnetic tape as described in black PDS manual when signing log book, put 'FORTH' in comments column. To close down see PDS manual. In the following documentation, the response OK will be assumed after all commands.

[2] Options
 [2.1] List of Options
 Select an option; loading is described under the paragraph number given. Subsequently loading another option removes the programs for the first option.*

| <u>Option name</u> | <u>Brief description</u> | <u>Go to</u> |
|--|--|--------------|
| Bigscan | Scan a plate; store data on tape, or on disk Maximum row size is 3000 pixels Additional option-result may be averaged down by factor of 2^n in x and y, while scanning. | [3] |
| Taylor | Scan a plate, store data on disk Maximum row size is 512 pixels Align a plate, by creating coordinate transform between it and a reference plate. Graphic output; range of additional sub-options are available, eg averaging down in x and y as in Bigscan, adding, subtracting and dividing two images; contour plots, moving arrays between disk and tape, defect removal, disk file graphic display | [4] |
| Smoothing | Averaging digitized scan stored on disk, without shrinking array (boxcar convolution) | [5] |
| 3D | 3D plot of a scan stored on disk or tape (several versions) | [6] |
| Photometry | Fitting 2D Gaussians to digitized scan (on disk or tape) | [7] |
| <u>Additions</u> | | [8.1] |
| Very long scans | More than 3000 pixels/row | [8.2] |
| Multiple scans | Reflecting arrays in x and y axes | [8.3] |
| Array flipping | | [8.4] |
| Tape file record-by-record graphic display | | [8.5] |
| Coarse file | | [8.6] |
| Density histogram | | |

*Except for options covered in [8.1], [8.2] and [8.5] which are to be loaded 'on top of' BIGSCAN and [8.3], to be loaded with either BIGSCAN or TAYLOR. This is made clear in the relevant documentation.

[3] BIGSCAN

[3.1] If you wish to use the capability to average down in x and y, as opposed to simple scan onto tape, go to [3.3]

Otherwise:

[3.2] Loading BIGSCAN

κ BIGSCAN LOAD

κ is the maximum row size in pixels that you will be using. It must be ≤ 3000 . For larger scan lines see [8.1] Wait for the OK before proceeding.

Go to [3.4]

[3.3] Loading BIGSCAN with Averaging Capability

When averaging, a scan of ρ rows of κ pixels each will be stored in an array on tape of ρ/a rows with κ/a pixels each, where $a = \frac{n}{2}$ is the value of a constant called AV. However, working storage must be reserved for the full κ pixels.

κ BIGSCAN LOAD

κ is the maximum row size in pixels, before averaging, that you will be using. It must be ≤ 2800
Loads averaging routines.

Continue

435 LOAD

[3.4]

Scanning

Set up over desired scan

Scan may be specified by center or any corner.

If you wish to zero PDS coordinates at this point,

ZILCH

If you wish to scan in x no special command is needed, this is default.
However if you scan in y and later wish to scan in x, you must then use

HORIZONTAL

If you wish to scan in y VERTICAL

HELLO ? JCM LOAD 26/07/78 33:52 OK
19:48 UT OK
1024 BIGSCAN LOAD OK
ZILCH OK
3 MM/S OK
0.000 ! OK
PLATE 64 64 WIDE 10 10 EXDY! OK
CENTER OK
BEHIND OK
NTP OK
WRITE JCM TEST SCAN 1970 JUL26 OK

Example of a simple scan using BIGSCAN

VERTICAL OK
5 VEL ! 1.000 ! OK
PLATE 32 20 WIDE 12 5 EXDY! OK
LL CORNER OK
WRITE VERTICAL SCAN, 000=1 OK

Example of a further scan using non-standard parameters

To set PDS velocity

Either
or

n MM/S
m VEL, !

Sets velocity in mm s⁻¹
Sets velocity in PDS units
(m = 255/60n)

To set scan mode:

Either

Ø ' ODD !

Default. Scans in l direction only
(increasing X, or Y if VERTICAL)
Scans back-and-forth, and flips rows in
working storage so final result is same

or

l ' ODD !

Then set up scan parameters as follows:

PLATE
ρ κ WIDE

makes PLATE the current image (see [0.8])
Size of scan; ρ rows with κ pixels in each
row.

δx δy DXDY!
Sets parameters $\tilde{WN}(=CN)$

Sets sampling interval in microns
dx is interval in x-coordinate
dy is interval in y-coordinate
(sets parameters DX)

STATUS

If you wish to check
the above parameters,

Ignore output column marked 'ORIGIN'

If you wish to define the scan by its center:

Position over center
CENTER

Stores origin (ie coordinates of most -ve
Corner) by setting parameters CQ_2 using
current PDS position (x,y) and WN and DX,
such that $x_0 = x - \frac{1}{2} w_x \delta x$; $y_0 = y - \frac{1}{2} w_y \delta y$
where (w_x, w_y) are (κ, ρ) or (ρ, κ) depending
on scan direction.

Go to [3.5]

Otherwise define the appropriate corner

Position over corner

Either
or

UR
UL

Upper right corner; default
Upper left corner

Lower right corner
Lower left corner
These define corner of scan, as seen on PDS viewing screen, not on plate itself. They modify the signs of PLATE parameters DX Stores current PDS coordinates in PLATE parameters CØ. This position will be the first pixel ($\kappa = \emptyset$; $\rho = \emptyset$) of the scan.

LR
LL
CORNER
Then

[3.5] Scanning

For scanning to disk, see the description under TAYLOR LOAD, section [4.7.2]

If tape is a new tape

REWIND
NTP

Rewinds tape, initializes variables, Writes 1 record and 2 filemarks, then backspaces twice. The record is for compatibility with RT tapes which have no filemark before the first file

Otherwise the tape may be positioned anywhere before the final double filemark.

Then: If averaging option is loaded and you wish to use it on this scan, go to [3.6]

Otherwise

WRITE (identifier)

- (identifier) is a string of ≤ 40 characters (including spaces). Follow this by carriage return.
- (i) Tape will run forward until a double filemark is found, and will then backspace so that the second filemark is overwritten.
 - (ii) The scan will be performed using the PLATE parameters set previously. The data will be written on to tape in Scansalot format.
 - (iii) After the scan is completed, a double filemark will be written at the end of the file, and the tape will then be backspaced twice, so that a second WRITE will not require any repositioning of the tap

(For an example of user commands for a scan using BIGSCAN, see p)

Go to [3.7]

[3.6]

Averaging to tape

To average on $\kappa \times \rho$ scan into an array of size INT ($\kappa/a \times \rho/a$)

a ' AV !
R-UP TURN-OFF

AVWRITE (identifier)

Stores a ($= 2^n$) in constant AV.
Prints out warning message on other terminal
(VDU or DEWRITER if either is switched on.
The other terminal cannot be used while
AVWRITE is in use since this uses the other
terminals' program area as working storage.
As WRITE (see [3.5]) but averages.
eg If there are 18 pixels/row and AV=4, 1st
4 rows of 18 pixels each will be averaged
down to one row of 4 pixels and then written
on to tape, similarly for subsequent sets of
4 rows, until the 16th row, $\kappa, \rho \geq 17$ are
ignored.

Continue

3.7]

Further scans

You may return to section [3.4]; remember to type PLATE so that any parameters modified will be those of PLATE.
Reset any parameters by using MM/s, WIDE, DXDY! VERTICAL, ODD. You may leave any of these unchanged. Then
specify origin using CENTER or UR, UL, LL, LR CORNER

3.8]

To Read Files Created by BIGSCAN

R-DOWN

R-DOWN lowers location of working storage
to program area and 2nd terminal, necessary
due to a hardware fault causing crashes
when reading to high memory. Before using
AVWRITE after this, use R-UP which raises
working storage again.

Go to [4.12] (Reading tape files)

[4]

[4.1]

Loading

TAYLOR LOAD

Loads PDS scanning routines, Tektronix graphics routines, least squares fit routines. STARS file routines, and disk file index. Max no of pixels/scan row is 512.

[4.2]

Sub-options

There are a number of additional facilities which may be loaded as well as the programs in TAYLOR LOAD. Only one of these options may be loaded at a time, subsequently loading another one will remove the first.

Either (i)

PROCESSING LOAD

Loads programs to aver in i,j and to add, subtract, divide and normalize two disk files

Go to [4.3]

or (ii) Either
or

REMOVAL LOAD) see [4.9] before
SLICING OPT LOAD) loading

Loads programs to display data from a file line by line, and to remove defects by linear interpolation.

Go to [4.9]

or (iii)

CONTOUR LOAD

Loads programs for contour plot of disk file.

or (iv)

SCANSALOT LOAD

Facility to move disk files to tape, and vice versa, read tape files etc. Also scans direct to tape

Go to [4.12]

4.3]

If you wish to align a plate before a scan, so that corresponding arrays may be obtained from different plates for the purpose of adding two plates of same object, or for the normalization of electronographs using a photocathode map go to [4.4]

Otherwise If you just wish to scan to disk without using the plate alignment facility, go to [4.7]

ALIGN

Description of method: 3 or more points (fiducial marks or stars) are entered into a file, called the STARS file, followed by the scan origin coordinates. A scan is done. Corresponding points are then entered in file for a second plate. Least-squares fitting is used to determine angle of rotation needed for the plates to be aligned (the translative component of the coordinate transformation is taken care of by the software but the rotation must be made zero by physically rotating the stage.) After stage is rotated, re-enter points and iterate until angle is acceptably small. (Usually 2 or 3 iterations are enough.) Then the corresponding scan is done.

(NOTE: the stars file is stored on disk so the data will not be lost until they are overwritten by the user.)

Continue

LINK-CLEAR

Clears out stars file.

If there is a set of data already in it (eg you are aligning different sets of plates) omit this instruction, and place the new data later in the file

Position over the first fiducial mark (or star)

If it is a star, you may wish to fit gaussians to it to determine the center point.

The word 'mark' will be used (for this section and the next) to mean fiducial mark or star.

(name) is ≤ 6 characters long and is for user identification only, not used in programs. This will so a scan 640 μ long (64 pixels) in Y (center of scan line is the current PDS position), fit gaussian to the scan line data, plot the data points and the fitted curve on the Tektronix screen. Do similar scan in X at Y value determined by first fit, fit gaussian, type out coordinates of fitted center and enter them and (name) into record 1 of STARS file. If you are keeping old data in records 1 through n of stars file, you can start at record n + 1 thus N+1 ENTER (name)

(name) is ≤ 6 characters long and is not used in programs. Enters (name) and current PDS coordinates


If gaussians to be fitted:

1 ENTER (name)

Otherwise

Center over mark
1 \$ENTER (name)

TAYLOR LOAD OK
PROCESSING LOAD OK
LINK-CLEAR OK
1 ENTER SMALL



Example of ENTER

ROYAL GREENWICH OBSERVATORY

| # | NAME | X | Y | STARS FILE | X | Y |
|---|--------|--------|--------|------------|--------|---|
| 1 | FLU | -24778 | 2 | -24767 | -11 | |
| 2 | FRO | -24634 | -24620 | -24627 | -24604 | |
| 3 | FRL | -9478 | -11933 | -9547 | -11918 | |
| 4 | GALAXY | 1426 | -4209 | 1356 | -4182 | |
| 5 | STAR1 | -22561 | -5237 | -22626 | -5310 | |
| 6 | STAR2 | -21128 | -10525 | -21180 | -18498 | |
| 7 | STAR3 | -10825 | -11999 | -10891 | -11983 | |
| 8 | HII#1 | | | | | |

OK

These columns not used

Example of STARS ALL

into record 1 of STARS file. If you are keeping other data in the first n records of STARS file, you can start at record n+1, thus n+1 $\$$ ENTER (name)

Then
If you wish to check contents of STARS file at any time
STARS ALL

Follow by two carriage returns. Page will automatically be cleared when return key is pressed for a second time.

(or n+2 ENTER (name) if you started at n+1)
(or n+2 ENTER (name) if you started at n+1)

Proceed to next mark
2 ENTER (name)
2 $\$$ ENTER (name)

Then carry on, using ENTER, or $\$$ ENTER to enter consecutive records of the STARS file, until all marks have been entered. Some marks may be entered using ENTER and others using $\$$ ENTER.
After all marks have been entered, move the PDS to the position of the scan (or the first scan, if more than one) that you are going to do, and enter the position in the STARS file in the next record, using $\$$ ENTER (eg 4 SENTER SCAN). Enter all the scans that you wish to re-scan on the second plate in consecutive records of the STARS file using $\$$ ENTER. The point that defines a scan can be its center or its center or any of its corners.

To do the scans go to [4.7] or [4.12]
for tape and then after scanning go to [4.6].

[4.6.]

Aligning second plate

You wish to do a scan or scans which will correspond pixel for pixel to a scan done on another plate, whose fiducial or stars have been entered in the STARS file. The new plate must be already roughly aligned by eye with the first one, within a few degrees, (by aligning fiducial marks or stars it should be possible to rotate the plate to within about 5° or the other: this is quite adequate.) The program uses the approximation $\sin \theta \approx \theta$.

[4.6.1]

If gaussians to be fitted 1 FIRST
Position over first mark ENTERed on first plate

Initializes rotation variables. Scans in Y and X as in ENTER, displays results as in ENTER, but stores final coordinates in second set of X, Y columns in STARS file. If for first mark on first plate you used a record other than record 1 (ie n+1 ENTER (name)) use the same record here (ie n+1 FIRST).

Otherwise 1 $\$$ FIRST

Initializes rotation variables Enters current coordinates in second set of

X,Y columns in STARS file.

The record number used should be the same as for the corresponding mark on the first plate, ie if (n+1) ENTER (name) or n+1 \$ENTER (name) was used, n+1 \$FIRST should be used here.

NOTE: The records in the file used for the alignment must be consecutive, but if the first mark cannot be found, find another and use m FIRST where m is the appropriate record number, then type \emptyset R## ! (or n R## ! if first mark is record n+1) before using NEXT (see below).

Then [4.6.2] NEXT

PDS will move to a position corresponding to next record in the STARS file, plus an offset of the coordinate origin calculated from the 2 sets of coordinates in the record for which FIRST or \$FIRST was used; ie where the next mark ought to be if the angle of rotation were zero.

Move the PDS to be correctly over the mark

In the first iteration the error may be quite large but in later iterations the position reached after NEXT should be close enough to require only a little 'tweaking up'

If gaussian fit required

OK

Scans, fits, displays, stores in STARS file (second set of columns). Record number was incremented by NEXT.

Otherwise

\$OK

Stores current coordinates.inSTARS file (second set of columns). Record number was incremented by NEXT, so will be correct automatically

Then loopback to [4.6.2] until all marks have been stored.

Then

n m ROTATE

STARS filerecords n through m inclusive are used to fit a transformation between the two sets of coordinates. If first record was record l and 3 marks were entered,
1 3 ROTATE

ROTATION

X2 Y2 ANGLE

23.4 -14.2 3.760

12.0 15.0 0.218

137

Rotate PDS stage through the appropriate angle

Output from ROTATE

The number immediately under ANGLE represents the angle of rotation in degrees X2, Y2 are offsets and second, third rows are residuals. The only number the user needs is the angle.

+ve angle \equiv rotate clockwise (plattten)

If angle is large, accuracy in rotation is not important due to the approximation in the transformation

Iterate by returning to [4.6.1], repeat until angle is sufficiently small.

Then

NEXT

This will take the stage to the appropriate scan position (ie the position corresponding to the scan location entered for the first plate). The second scan may now be done into a different disk or tape file. (see [4.7] or [4.12]). If several scan positions were stored in STARS file for the first plate, repeating NEXT will move stage to the corresponding positions, so all the scans may be done in turn.

[4.6.3]

Additional ALIGN commands

y x GO

n 1GO

n 2GO

Where y, x are 32-bit PDS coordinates; moves stage to these coordinates.

moves stage to coordinates given by record n of STARS file (1st set of coordinates).

moves stage to coordinates given by record n of STARES file (2nd set of coordinates).

Graphics mode: The data points for the gaussian fit in ENTER, FIRST, OK are plotted out as crosses. ("POINTS mode") Alternatives are listed below, using the word sets the plotting mode until another of the words is used or the graphics are re-loaded.

POINTS

FOLLOW

HISTOGRAM

The commands AXES and BIG GRID will plot axes and a grid on the screen.

[4.7] Scanning

[4.7.1]

Position PDS over point defining first scan (if ALIGN has been used, this will already have been done.) This point may be the center or any corner.

If you wish to zero the PDS coordinates at any time, this can be done by the command ZILCH or by pressing the appropriate buttons on the console displaying the coordinates.

If you wish to scan in Y

VERTICAL

Sets Y-scan flag

To return to scanning in X later,

HORIZONTAL

Clears Y-scan flag

Otherwise: Scanning in X is default; it is not necessary to type HORIZONTAL unless VERTICAL has been used earlier.

To set PDS velocity

Either
or

n MM/S
m VEL !

Sets velocity in mm s⁻¹
Sets velocity in PDS units (m = 255/60 n)

To set scan mode:

Either
or

0 ' ODD !
1 ' ODD !

Default. Scans in 1 direction only (increasing X, or Y if VERTICAL)
Scans back-and forth, and flips rows in working storage so final result is same.

Then set up scan parameters as follows:

PLATE
ρ κ WIDE

Makes PLATE the current image (see [0.8]).
Size of scan; ρ rows with κ pixels in each row. Sets parameters WN, and sets CN equal to WN.

δx δy DXDY!
Sets sampling interval in microns δx is interval in x-coordinate, δy is interval in y-coordinate (sets parameters DX)

If you wish to check the above parameters

STATUS

If you wish to define the scan by its center

Position over center
CENTER

[4.7.2]

Otherwise define the appropriate corner.

Position over corner

Either
or UR
or UL
or LR
or LL

CORNER

Then

If you are scanning to tape using SCANSALOT LOAD, return to [4.12]

Otherwise go to [4.7.2]

Scanning to disk

[4.7.2]

SOURCE

Ignore output column marked 'ORIGIN'

Stores origin (ie coordinates of most -ve corner) by setting parameters $C\emptyset$, using current PDS position (x,y) and \widetilde{WN} and \widetilde{DX} , Such that $x_o = x - \frac{1}{2} w_x \delta x$

$$y_o = y - \frac{1}{2} w_y \delta y$$

where (w_x, w_y) are (κ,ρ) or (ρ,κ) depending on scan direction.

Upper right corner: default

Upper left corner

Lower left corner

Lower left corner

These define corner of scan, as seen on PDS viewing screen, not on plate itself. They modify the signs of PLATE parameters DX .

Stores current PDS coordinates in PLATE parameters $C\emptyset$. This position will be the first pixel ($i = \emptyset, j = \emptyset$) of the scan

Specifies PLATE (the current image) as the source for the pass (scan operation). The parameters set up in PLATE will be copied to INPUT. SOURCE leaves INPUT as the current

ROYAL GREENWICH OBSERVATORY

IMAGES INDEX

| # NAME | LOC'N | ORIGIN: X | Y | DIR. | CORN. | SIZE: X | Y | +X | +Y |
|-----------|-------|-----------|--------|------|-------|---------|-----|----|----|
| 1 CART2 | 485 | 0. | 0. | H | UR | 128 | 128 | 10 | 10 |
| | 516 | 640. | 640. | | | | | | |
| 2 CART3 | 520 | -1920. | -1920. | H | UR | 64 | 64 | 60 | 60 |
| | 527 | 0. | 0. | | | | | | |
| 3 COPY | 470 | 0. | 0. | V | LL | 64 | 64 | 10 | 10 |
| | 477 | -320. | -320. | | | | | | |
| 4 JCM | 400 | -2560. | -2560. | H | UR | 128 | 128 | 40 | 40 |
| | 431 | -2 | 0. | | | | | | |
| 5 LONG | 600 | 0. | 0. | V | UR | 512 | 4 | 10 | 10 |
| | 603 | 20. | 2560. | | | | | | |
| 6 M0001 | 457 | -1000. | 1000. | H | UL | 32 | 32 | 10 | 10 |
| | 458 | -1160. | 1160. | | | | | | |
| 7 N1068#1 | 430 | -430. | -501. | H | UR | 50 | 50 | 20 | 20 |
| | 434 | 14. | -1. | | | | | | |
| 8 PK | 461 | 0. | 0. | H | LR | 64 | 64 | 10 | 10 |
| | 468 | 320. | -320. | | | | | | |
| 9 TEST | 459 | -1000. | 1000. | H | UL | 32 | 32 | 10 | 10 |
| | 460 | -1160. | 1160. | | | | | | |
| 10 TWIN | 432 | 0. | 0. | V | UR | 64 | 64 | 10 | 10 |
| | 439 | 320. | 320. | | | | | | |

OK

Which corner?
(COUNTER=UR)

X pixels Y pixels
(K) (e)

Horizontal or
Vertical scan?

X origin Y origin

X center Y center

Starting block
Finishing block

(Name)

dy
microns

Example of LISTFILE

image (ie use of WIDE etc will modify INPUT parameters, not those of PLATE.

LISTFILE (cr) (cr)

LISTFILE must be followed by two carriage returns, it automatically clears the page (on the Tektronix) It prints out index of data for current disk files (see example attached) If the index takes more than 1 page, press any key except PAGE to clear page.

This is highest number in column marked LOC'N

Make an attempt to fill up empty disk areas.
1 Block takes 512 pixels (32 x 32 array = 2 Blocks.)
n is starting block number ($500 \leq n \leq 2435$)
(filename) is a string of ≤ 8 characters.
This creates a file accessed by the command
FIND (filename)

If there are old files in the index that are no longer used,

REMOVE (filename)

Removes entry from index; Files should be removed as soon as they are finished with. This command can be used at any time.

the averaging capability

If you have done a PROCESSING LOAD (or do one at this point) and wish to use the
for this scan, go to [4.7.3]

Otherwise

FIND (filename) EQUAL DESTINATION

FIND (filename) selects the file (the one just created, usually) EQUAL copies INPUT's parameters to DISK and DESTINATION specifies DISK as the OUTPUT image, copying parameters form DISK to OUTPUT.

PASS

Passes INPUT to OUTPUT. Since INPUT is PLATE and OUTPUT is DISK, PASS will do a scan according to the parameters set up, and store the data on disk, Remember to select AUTO on PDS

If you are using ALIGN and now wish to align second plate, go to [4.6].

If you wish to do further scans, go to [4.7]

For further capabilities in PROCESSING LOAD see [4.8]

[4.7.3]

Averaging

```
FIND (filename)
a ' AV !
AVBY DESTINATION
```

Sects disk file
Stores a (number to be averaged by,
 $a = 2^n$) in AV
AVBY takes INPUT parameters, multiplies or
divides by AV where appropriate, and puts them
in DISK's parameters. DESTINATION specified
DISK as the OUTPUT (copies DISK's parameters
to OUTPUT)

AVERAGE

Modified PASS - averages in i and j scans
a line, puts this in working storage. Takes
first AV numbers, adds each into first loca-
tion of "sum buffer"; next AV numbers into
second location, etc. Before being added
in, numbers are divided by AV^2 2nd, 3rd...
AVth scan line added in and then sum buffer
stores as 1st row of output array, and sum
buffer zeroed. Repeated for every group of
AV scan lines, so final array will have
each number being the average of AV^2 numbers;
eg a scan 64 64 WIDE with 4 in AV will result
in an output array of size 16 by 16 with
sample size multiplied by 4 in each direction
(ie 10 10 DXDY! Will have an output dx, dy
of 40, 40).

If you are using ALIGN and now wish to align second plate, go to [4.6]

If you wish to do further scans, go to [4.7]

For further capabilities in PROCESSING LOAD see [4.8]

Further capabilities in PROCESSING LOAD - combination of disk images

4.8]

For two disk files (file 1) and (file 2) to be arithmetically combined to form a third file (file 3) use
the following instructions

```
FIND (file 1) SOURCE
FIND (file 3) DESTINATION
```

(file 1) is identified with INPUT
(file 3) is identified with OUTPUT

(file 2) is identified with OPERAND by using one of the words.
 PASS combines data in INPUT with data in OPERAND and places the result in (file 3)

```

( ADD
( SUBTRACT
( DIVIDE
PASS
  
```

The commands defining the OPERAND are

```

ADD giving (file 3) = (file 1) + (file 2)
SUBTRACT giving (file 3) = (file 1) - (file 2)
DIVIDE giving (file 3) = NUM x (file 1)/(file 3)
  
```

where NUM is a variable (default 1000) which may be set by n NUM :

Each element in the INPUT array is combined with the corresponding element in the OPERAND array row by row in working storage. The INPUT and OPERAND files are unaltered.

If you wish to overwrite (file 1) with the result then

```

FIND (file 1) SOURCE DESTINATION
) ADD
FIND (file 2)) SUBTRACT PASS
) DIVIDE
  
```

If, however you wish to overwrite (file 2) with the result, then

```

FIND (file 1) SOURCE
FIND (file 2) DESTINATION
ADD )
SUBTRACT) PASS
DIVIDE )
  
```

[4.8.2]

Sum and average of an array

SBAR will calculate the average value of an array and both print it out and store it in the variable NUM. The sum of the array is stored in the double precision (32-bit) variable XSUM, to see value: XSUM 2@ D. To subtract a constant background value n throughout the calculation (while not altering the array itself) use the command

```
n BKG !
```

before using SBAR; default BKG is zero, after changing it, it will stay at current value until altered again or until PROCESSING is reloaded.

```
FIND (filename) SBAR AVERAGE VALUE = 1025
```

Calculates sum and average

[4.8.3]

Normalization of arrays

For example, for removing variations due to photocathode sensitivity on an electronograph, divide scan array by scan photocathode map and multiply through by average of latter scan.

Suppose first scan is in a disk file called STAR, the photocathode map in a file called SKY, and you wish to place the result in a file called FINAL, which you first need to create with size equal to those of STAR and SKY, starting at block n.

```
LISTFILE
FIND STAR
      n INTER FINAL
```

See [4.7.2] for details
Makes disk current image with STAR's parameters.
Enters FINAL in index file with parameters equal to STAR.

```
FIND SKY SBAR AVERAGE VALUE = 952
```

Stores average value of SKY by NUM.

```
FIND STAR SOURCE FIND FINAL DESTINATION
```

Divides STAR by SKY

```
FIND SKY DIVIDE PASS
```

multiplying through by NUM and putting the result in FINAL.

If you wish to do further DIVIDE's with NUM at the default value

```
1000 NUM !
```

[4.8.4]

NOTES (i) If you wish to use only part of the image, ie a subset of window, see notes on windowing in [4.11]

(ii) Words like ADD, DIVIDE, SUBTRACT do not do anything themselves except assign an operation called PROCESS which is executed for each row in PASS. At the end of a PASS, a default operation which does nothing is re-assigned so that subsequent PASSES will not perform any operation. This default operation is called NOTHING; using it immediately after ADD, etc will cancel the effect of the latter.

These words also assign the current image to be the OPERAND (they include the word SECOND).

[4.9]

REMOVAL LOAD/SLICING LOAD

SLICING OPT LOAD loads basic display and defect removal routines. (disk files only)
REMOVAL LOAD loads the above, plus routines to

- (1) Makes PDS move to appropriate location for direct visual examination of defects
- (2) Make available a file of defect coordinates
- (3) Fit gaussians to scan lines

If you only want the basic SLICING routines

SLICING OPT LOAD
Go to [4.9.1]

Loads limited
options

Otherwise, for REMOVAL

REMOVAL LOAD
Go to [4.9.1], continuing through to [4.9.2]

Loads complete
package

[4.9.1]

SLICING

To view an array line by line:-

FIND (filename)
FIRST
N



Select appropriate disk file.
Displays first line, sets graphics scalings and row sets no R to \emptyset .
Displays next line. Repeat N for consecutive lines clearing pages as required.

[If response is OK with no lines plotted, file is complete]

Subsidiary Display words:-

PLOT
n ' R + !
(n-1) ' R : N
m n X SIZE
m n Y SIZE

Replots current line. Useful after clearing page.
Skips n lines
Plots line n
Resets x-scale so that pixels $m \leq i \leq n$ are plotted
Resets Y-scale so that density values m through n are plotted. (Default to 0.4000)

Sets Plotting mode + + + + +
 Sets Plotting mode 
 Sets Plotting mode 
 (Default.)

Puts up cursors

Program does linear interpolation between 2 points.

POINTS

HISTOGRAM

FOLLOW

To removed defect, plot line and

FIX

Position cursor on L.H.S. of defect

Tap any key (except RETURN)

Position cursor on R.H.S. of defect

Tap key

[New interpolated line is stored in the disk file and replotted on screen.
 Old line is lost].

These are all the available words in SLICING.
 If you have loaded REMOVAL continue to [4.9.2]

[4.9.2] REMOVAL

This assumes that source of disk file is still on PDS and co-ordinate system has not been changed.

[4.9.2.1]

Plot lines until possible defect appears

If desired clear page and PLOT

If you wish to enter removed-defects co-ordinate into DEFECTS file to initialize that file then

?DEFECTS (or) (or)

Displays cursor and contents of DEFECTS

LINK-CLEAR

Initializes file

Otherwise

NOX

Centre cursor over apparent defect

Puts up cursor

This position is not used in any subsequent defect removal and hence position of cursor is not critical.

Tap any key (except RETURN)

PDS will move to co-ordinate corresponding to that point and the user may decide from visual appearance of area whether removal of feature is desired.

If defect to be removed go to [4.9.2.2]
Otherwise go back to [4.9.2.1] and locate another defect candidate

[4.9.2.2] If defect to be entered into DEFECTS file go to [4.9.2.3]
Otherwise use FIX as described before in [4.9.1.1]

[4.9.2.3] NIX

Identical to FIX except, stores current PDS co-ordinates in DEFECTS file. (Do not move PDS between NOX and NIX otherwise wrong co-ordinates will be entered.)

To view DEFECTS file

?DEFECTS

Defects are entered sequentially numbered after last defect in file.

To remove further defects go to [4.9.2.1]

If a second plate is then set up using ALIGN so that the co-ordinate translation has been calculated using ROTATE

n @DEFECT

This will move PDS to position corresponding to coordinates in record n of defects file, with the translation added.

In this way the user can see if the same defects are present on different plates.

If no ROTATE has been performed since loading TAYLOR then this instruction will take PDS to nth defect of existing plate.

[4.9.2.4]

Fitting a Gaussian to a Scan Line

Plot Line

by using (N or PLOT)
APPROX FIT

Calculates 1st approximation
and does least squares fit

Plots fitted curve

Outputs values of parameters XØ

(value in column CENTER) is distance
from center of curve in units of
 $\frac{1}{10} \mu$ (i.e. if result is x, pixel number
in line is $\frac{x}{10} \delta x + \frac{1}{2}K$ where K is pixel
length of row)

| CURVE RESULT | CENTERING | | | |
|-----------------|-----------|------|------|------|
| | CENTER | FEAK | BASE | HW |
| | 1.52 | 3000 | 257 | 3.42 |
| | 0.00 | 0 | 0 | 0.00 |
| | 0 | | | |

PEAK Is height of the peak above the base
BASE Is the base. Units for these are the
same as the scan data

(Ø ≤ density ≤ 4095)

HW Halfwidth, units $\frac{1}{10} \mu$. Second line gives
residuals

Repeat FIT, CURVE, RESULT as required, for
example:

FIT CURVE FIT FIT CURVE FIT RESULT CURVE etc

Selects disk image to be plotted

[4.10] CONTOUR LOAD

[4.10.1]

FIND (filename)

If only a part (window) of the file is to be plotted, see [4.11] to do this.

SCREEN n WIDTH

sets up the screen to be n wide
by $\frac{3}{4}$ n high (in pixels)

So, for an m x m image -
 $\frac{4}{3}$ WIDTH Will set the scale so
that all the file will be plotted,
and the right hand quarter of the
screen will be blank.

m WIDTH would fill the whole screen but
only plot the first $\frac{3}{4}$ m rows

n SPACE :

Optional: default 1000
Sets contour interval to be n
(i.e. difference of n in the numbers
in the array.)

low high CLIP

Optional: default -32768 32767 CLIP
"low" and "high" are numbers defining
cutoff points, for example 200 4000
CLIP will treat all numbers <200 as
200, and all >4000 as 4000

DISK PLOT

Plots contour map

[4.10.2]

Contour plotting directly from tape


512 MAXIMUM SCANSALOT LOAD CONTOUR. LOAD

This unloads TAYLOR so if you
subsequently wish to use programs in
TAYLOR you must repeat TAYLOR LOAD.

n SKIPS TWIDE

Select tape file; described in [4.12]

Q 1200 CLIP 100 SPACE 1 DISK PLOT



OK


Example of a PLOT using CONTOUR

If windowing is desired, see [4.11]

SCREEN n WIDTH Described in [4.10.1]
 n SPACE ! low high CLIP " " "
 TAPE PLOT Plots contour map.
 PIXELS 0 1 2 . . . M . . . K-1

[4.11] Disk file manipulation

[4.11.1] Windowing



Given a 2D array A_{ij} of pixel dimension $\rho \times n$ (such that $0 \leq i < \kappa$ and $0 \leq j < \rho$) a subset B_{ij} of pixel dimension $a \times b$ (such that $m \leq i < (m+a)$ and $n \leq j < (n+b)$) may be accessed, provided m and $n \geq 0$, $(m+a) < \kappa$ and $(n+b) < \rho$.

This is used for plotting part of a large array, combining part of an array with a smaller array fitting a function to part of an array, etc.

Windowing words are:

n m ORIG

$\vec{W0}$: defines window origin
 at row n, pixel m

| | | | | | | | | | | |
|------|------|------|------|------|------|------|------|------|------|-----|
| 59 | 59 | 60 | 61 | 62 | 63 | 64 | 65 | 66 | 67 | 68 |
| 541 | 554 | 554 | 554 | 555 | 554 | 551 | 541 | 551 | 564 | 540 |
| 554 | 557 | 557 | 562 | 562 | 569 | 570 | 575 | 568 | 565 | 559 |
| 557 | 576 | 577 | 591 | 591 | 592 | 586 | 587 | 588 | 580 | 584 |
| 560 | 592 | 595 | 597 | 597 | 635 | 611 | 609 | 607 | 602 | 592 |
| 595 | 620 | 642 | 642 | 655 | 655 | 655 | 654 | 647 | 632 | 620 |
| 624 | 656 | 654 | 691 | 718 | 718 | 717 | 717 | 695 | 676 | 665 |
| 672 | 700 | 726 | 742 | 742 | 755 | 775 | 776 | 764 | 738 | 701 |
| 701 | 755 | 789 | 823 | 823 | 845 | 836 | 839 | 830 | 794 | 744 |
| 793 | 818 | 860 | 917 | 917 | 946 | 945 | 933 | 908 | 860 | 797 |
| 853 | 884 | 959 | 1030 | 1030 | 1077 | 1100 | 1047 | 995 | 937 | 858 |
| 880 | 976 | 1077 | 1177 | 1177 | 1245 | 1255 | 1204 | 1104 | 1015 | 914 |
| 943 | 1062 | 1165 | 1318 | 1318 | 1430 | 1440 | 1363 | 1320 | 1082 | 960 |
| 983 | 1116 | 1233 | 1470 | 1470 | 1591 | 1581 | 1475 | 1300 | 1133 | 988 |
| 1002 | 1152 | 1346 | 1532 | 1532 | 1685 | 1669 | 1515 | 1326 | 1154 | 996 |
| 1044 | 1162 | 1373 | 1600 | 1600 | 1849 | 1614 | 1473 | 1295 | 1118 | 980 |
| 1066 | 1183 | 1406 | 1827 | 1827 | 1956 | 1450 | 1356 | 1211 | 1064 | 939 |
| 1093 | 1227 | 1242 | 1948 | 1948 | 1276 | 1276 | 1203 | 1093 | 979 | 877 |
| 1115 | 1255 | 1643 | 2035 | 2035 | 1699 | 1699 | 1053 | 996 | 914 | 820 |
| 1130 | 1266 | 1719 | 2110 | 2110 | 1962 | 1962 | 927 | 894 | 831 | 760 |
| 1145 | 1291 | 1853 | 2233 | 2233 | 2059 | 2059 | 830 | 795 | 759 | 710 |
| 1160 | 1316 | 1929 | 2306 | 2306 | 2084 | 2084 | 755 | 729 | 702 | 671 |
| 1175 | 1341 | 2066 | 2384 | 2384 | 226 | 226 | 699 | 681 | 660 | 644 |
| 1190 | 1367 | 2166 | 2467 | 2466 | 676 | 660 | 651 | 644 | 629 | 612 |
| 1205 | 1391 | 2211 | 2541 | 2541 | 633 | 633 | 616 | 611 | 596 | 582 |
| 1220 | 1417 | 2311 | 2600 | 2600 | 594 | 585 | 587 | 579 | 570 | 549 |
| 1235 | 1443 | 2399 | 2667 | 2667 | 567 | 560 | 563 | 565 | 557 | 547 |
| 1250 | 1469 | 2490 | 2740 | 2740 | 549 | 540 | 542 | 549 | 546 | 536 |
| 1265 | 1495 | 2541 | 2811 | 2811 | 540 | 528 | 527 | 527 | 532 | 532 |
| 1280 | 1516 | 2631 | 2888 | 2888 | 540 | 518 | 516 | 517 | 516 | 519 |
| 1295 | 1541 | 2731 | 2966 | 2966 | 531 | 513 | 513 | 511 | 509 | 505 |
| 1310 | 1566 | 2831 | 3044 | 3044 | 531 | 513 | 513 | 511 | 509 | |

Example of VOILA

OK

Sets \vec{WN} , defines size of window to be b rows of a pixels each.
Sets $\vec{W0}$ to (0.0) and $\vec{WN} = \vec{CN}$ ie; makes the window the default value, the whole array.

[4.11.2] Digital display

FIND (filename) m n VOILA

Clears page and displays window of file with origin at pixel n, row m. Resets $\vec{W0}$, \vec{WN} to default at end.

[4.11.3] Examples of Copying Files

FIND (file 1) SOURCE
FIND (file 2) DESTINATION PASS
FIND (file 1) 32 32 WW SOURCE
FIND (file 2) 32 32 WW DESTINATION PASS

Copies (file 1) to (file 2)
Copies "top left" (1st 32 pixels of 1st 32 rows of file 1) to corresponding area of (file 2)

[4.12] SCANSALOT LOAD

Taping by SCANSALOT LOAD, uses the second terminal's working storage area so the second terminal may not be used while you are using the tape.

[4.12.1] To set up parameters for a scan to tape, go to [4.7.1]

If tape is a new tape (no files on it to be kept)

REWIND

Rewinds tape
Initialises tape variables

NTP

Writes 1 dummy record and 2 filemarks then backspaces twice. The record is for compatibility with RTSCAN tapes which have no filemark before the first file.

Otherwise ensure tape is positioned at some point before the final double filemark. Program will look for a double filemark before starting to write.

Then
WRITE (identifier) (cr)

(Identifier) is a string of ≤ 40 characters (including spaces). Tape will run forward until a double filemark is found, and will then backspace so that the second filemark will be overwritten. The scan will be performed using the PLATE parameters set previously. The data will be written onto tape in Scansalot format.

After the scan is completed, a double filemark will be written at the end of the file, and the tape will then be backspaced twice, so that a second WRITE will not require any repositioning of the tape.

The rest of this section [4.12] can be used without doing a TAYLOR LOAD, by doing a 512 MAXIMUM SCANSALOT LOAD instead.

[4.12.2]

Disk-to-tape pass

FIND (filename)
WRITE (identifier)

Selects appropriate disk file.
Looks for end of tape (double filemark) backspaces once and copies disk file onto tape. Finishes by writing a new double filemark and backspaces twice.

(Identifier) is ≤ 40 character string including spaces terminated by carriage return.

REWIND

Rewinds tape, initialises tape variables EOFS and REC.

n SKIPS

n > 0

Passes forward over n filemarks. Remains positioned immediately after the last (i.e. at the beginning of a file.)

n < 0

Goes backward over n filemarks then forward l space, so positioning tape at the beginning of a file

-1 SKIPS

Returns to beginning of current file.

TWIDE

Returns to beginning of current file, goes through file checking number of records and coordinates in header, uses these to set TAPE parameters \overline{WN} , \overline{CN} , $\overline{C0}$, \overline{DX} . Types first identifier, returns to beginning of the file.

T-INDEX

n SKIPS TWIDE is the taping equivalent of FIND (filename) for disk files.

Prints directory of files on tape; rewinds tape and goes forward until double filemark found, then backspaces twice.

n FILE

Positions at end of first record of file n, n is absolute file number with dummy header file being zero. However in some circumstances when interrogating individual records in a file, one goes past the end of file. Subsequent use of FILE may give the wrong result.

T-INDEX

- 0 DUMMY FILE ; BEGINNING OF TAPE
- 1 JCM TEST SCAN 1978JUL26
- 2 VERTICAL SCAN, ODD=1
- 3 ANOTHER VERTICAL SCAN
- 4 SIMPLE SCAN
- 5 AWRITE TEST

END OF FILES OK

| | | | | | |
|----|----|-----|----------|--------|---|
| 2 | 1 | 0 | 0 | 0 | 0 |
| 64 | 64 | 10 | 10-26860 | -3398 | |
| 20 | 32 | -12 | -5-14010 | -11834 | |
| 20 | 32 | -12 | -5-14010 | -11834 | |
| 64 | 64 | 10 | 10-26860 | -3398 | |
| 64 | 64 | 10 | 10 | 0 | |

WIN WNTT OX OXTT CD CD DT

Example of T-INDEX

LABEL

Reads and prints next records identifier.

n RECS

Moves forward n-1 records and reads next record

1 RECS reads next record

2 RECS reads next record but 1

0 RECS reads record just read

-1 RECS reads record before one just read

SHOW

Types identifier and dumps scan data currently in working storage. To be used after RECS.

[4.12.4] Tape to disk pass

n SKIPS TWIDE

Selects file, prints identifier, sets TAPE parameters

SOURCE

Copies TAPE parameters to INPUT

LISTFILE

See [4.7.2] for explanation, lists disk file

n INTER (filename)

See [4.7.2] for explanation, enters disk fi

FIND (filename)

Selects disk file

DESTINATION

Copies DISK parameters to OUTPUT

FIND (filename)

Selects disk file

EQUAL DESTINATION

Sets DISK and OUTPUT parameters

PASS

Copies tape file to disk file

Otherwise

Then

[5]

SMOOTHING (boxcar convolution)

[5.1]

Averaging in i and j by constant factor a ; disk arrays only.
 Description: let contents of pixel i , row j in array A be A_{ij}

A is source array, B is smoothed array
 Each array consists of ρ rows of K pixels each
 $[x]$ means integer part of x

Action of program:

- $B_{ij} = A_{ij}$ for (i) $\emptyset \leq j < [\frac{1}{2}a]$, $\emptyset \leq i < K$
- (ii) $\emptyset \leq i < [\frac{1}{2}a]$, $\emptyset \leq j < \rho$
- (iii) $K - [\frac{a+1}{2}] + 1 \leq i < K$, $\emptyset \leq j < \rho$

Partially smoothed s.t. $B_{ij} = \frac{1}{a} \sum_{l=j-[\frac{a}{2}]}^{j+[\frac{a}{2}]-1} \sum_{k=i-[\frac{a}{2}]}^{i+[\frac{a+1}{2}]-1} A_{kl}$

for $[\frac{1}{2}a] \leq i < K - [\frac{a+1}{2}] + 1$, $\rho - [\frac{a+1}{2}] + 1 \leq j < \rho$

Smoothed s.t. $B_{ij} = \frac{1}{a^2} \sum_{r=j-[\frac{a}{2}]}^{j+[\frac{a+1}{2}]-1} \sum_{s=i-[\frac{a}{2}]}^{i+[\frac{a+1}{2}]-1} A_{ksl}$

for $[\frac{1}{2}a] \leq i < K - [\frac{a+1}{2}] + 1$, $[\frac{1}{2}a] \leq j < \rho - [\frac{a+1}{2}] + 1$

Inside fully smoothed region:

For odd values of a -

B_{ij} is the average of the $a \times a$ square of points in array A with center (i, j)

For even values of a -

B_{ij} is the average of the $a \times a$ square of points in array A with center $(i - \frac{1}{2}, j - \frac{1}{2})$ so for even a the image is translated by $\frac{1}{2}$ pixel in i and j . If an array is to be repeatedly smoothed by an even factor it is recommended that the array be 'flipped' each time before smoothing so that $A_{ij} \rightarrow A(k-i-1) (\rho-j-1)$. Thus the successive translations will cancel instead of adding. How to achieve this is described below.

[5.2]

Tape files

This program only works on disk files, to smooth a tape file first copy it to disk, either 512 MAXIMUM SCANSALOT LOAD or TAYLOR LOAD SCANSALOT LOAD, and then go to [4.12.4] for instructions on tape-to-disk pass.

[5.3]

Loading

a SMOOTHING LOAD 512

a is the full width of the boxcar in p and k ($a \geq 2$).
Output is maximum row size allowed.
This is 512 or less depending on value of a. For $a > 8$ the buffer sizes must be < 512 to fit in **ccre**.

[5.4]

Disk files

Allocate all disk files you will require using

LISTFILE

FIND (file 1)

Lists images index, for explanation see [4.7.2]

Sets current image parameters to those of (file 1). (file 1) is the source file i.e. the one you wish to smooth.

n INTER (filename)

n is block number, select as described in [4.7.2]. INTER allocates a file called (filename) the same size as (file 1). Use INTER to allocate all the files you will need.

[5.5] Smoothing

(i) Single smoothing

FIND (file 1) SOURCE

Defines (file 1) as the INPUT file, i.e.: file to be smoothed

FIND (file 2) CONVOLUTE

Copies (file 1) to (file 2) and smooths (file 2). CONVOLUTE is equivalent to DESTINATION PASS \$CONVOLVE (see below).

(ii) Smoothing into same file

FIND (file 1) \$CONVOLVE

Smooths (file 1), overwriting original data. If you wish to save original data, copy it into another file or use CONVOLUTE as in (i) and then FIND (file 2) \$CONVOLVE.

FIND (file 1) \$CONVOLVE \$CONVOLVE \$CONVOLVE

Repeatedly smooths (file 1).

(iii) Flipping and then smoothing

This is for use with repeated smoothing by even values etc.

FIND (file 1) SOURCE

(file 1) will be overwritten later, to save original first copy it.

FIND (file 2) RCONVOLUTE

Copies (file 1) to (file 2) reflecting it on the way so that pixel (i,j) of (file 1) is put in pixel (k-i-1,p-j-1) of (file 2), i.e. reflected in the center point of the array.

Then smooths (file 2).

Flips (file 2) into (file 1) and smooths (file 1).

SOURCE FIND (file 1) RCONVOLUTE

SOURCE FIND (file 2) RCONVOLUTE
SOURCE FIND (file 1) RCONVOLUTE
SOURCE FIND (file 2) RCONVOLUTE

Repeatedly smooths array, with the
 $\frac{1}{2}$ -pixel translation caused by each
smoothing cancelling.

If the final file is (file 1), the array will be the right way round, but
if it is (file 2), one more flip must be done to correct it.

Takes current file, flips into (file 1)

SOURCE FIND (file 1) RC-FLIP PASS

6] 3D plot

There are 3 variants loaded by:


Either 3D LOAD
or 3D/J LOAD
or 3DSLW LCAD

Then if tape files are to be plotted: SCANSALOT LOAD


Basic FORTH 3D plot
Hidden line 3D plot
Alternative hidden line plot
Loads taping routines

6.1] Form of plot

Orientation of axes



X represents pixels/row
Y represents density units (numbers in array)
Z is direction of decreasing row number i.e.
later rows are plotted towards the more
negative Z.



n m X SIZE sets the X scale, n being the abscissa of the left hand side of the screen, m being the abscissa of the right hand side.

Thus \emptyset WN @ X SIZE will cause a plotted row to fill the screen from left to right. (WN @ gives the value of K)

n m Y SIZE sets the Y scale

n m Z SIZE controls the spacing between the rows.

Default values: \emptyset 32 X SIZE

\emptyset 8192 Y SIZE

- 64 64 Z SIZE

The pitch and yaw of the image can also be altered.

Pitch inclines the Z axis from the horizontal

Yaw inclines the Y axis from the vertical, in the plane normal to Z.
Default pitch is 60° , yaw is zero.

[6.2]

Commands

Commands for the three variants are identical except in the case of PITCH and YAW.

SCREEN ρ K WIDE

No more than ρ rows and no more than K pixels/row will be plotted.

For a disk file: FIND (filename)
For a tape file: n SKIPS TWIDE

Then [6.2.1]

If only a part of the file is to be plotted:

n m ~~ORIGIN~~
b a WW
n m X SIZE
n m Y SIZE
n m Z SIZE

If 3D variant being used

S c PITCH
S c YAW

Otherwise (ie 3D/J or 3DSLOW variants)

nn . nn PITCH
nn . nn YAW

Selects disk file to be plotted.

n is relative file number;
see [4.12] for explanation

see [4.11] for explanation

see [4.11] for explanation

Sets up X axis co-ordinates
n = abscissa of L.H.S.
m = abscissa of R.H.S.
(Default 0 32 X SIZE)

(Default 0 8192 Y SIZE)


(Default -64 64 Z SIZE)
For explanation see [6.1]

Default 866 500 PITCH


Default 0 1000 YAW
S = 1000 sin θ c = 1000 cos θ

Default 60.00 PITCH

Default 00.00 YAW
nn.nn is angle in degrees. There
must be exactly 2 decimal places -
do NOT omit trailing zeros.



a PLOT using 3DSLOW



DK

Q Plot using SP/5

Then For a disk file: DISK PLOT
 For a tape file: TAPE PLOT

Then to replot using different window, scales, pitch or yaw: go to [6.2.1]

[7] Photometry

Description: Least squares fitting of a 2-D gaussian to scan data from a disk file. The file may be windowed (see [4.11]).

Consider a 2-D array of data A_{ij} where $0 \leq i < k$ and $0 \leq j < p$.

Fitting Parameters

| | | | | | |
|--------------|----------------|---|----------------------------|---|---|
| $X\emptyset$ | $= 6400 \ i/k$ |) | $(X\emptyset, Y\emptyset)$ |) | is centre of gaussian |
| $Y\emptyset$ | $= 6400 \ j/p$ |) | δ_j / δ_i |) | δ_j / δ_i = ratio of sample intervals in microns of the j and i directions. |

Parameter definition

PEAK Height of peak above base (units as in A_{ij})
 BASE Base (units as in A_{ij})
 HX Half width at half height in i direction (units as in $X\emptyset$)
 HY Half width at half height in i direction (units as in $Y\emptyset$)
 SR Correlation. Stored as integer. Value multiplied by 1000.

Function:

$$f(i,j) = (BASE) + (PEAK) 2^{\xi} \leq C$$

$$\text{where } \xi = \left(\frac{X-X\emptyset}{HX} \right)^2 + \frac{SR(X-X\emptyset)(Y-Y\emptyset)}{HX.HY} + \left(\frac{Y-Y\emptyset}{HY} \right)^2$$

Clears page, plots file
 Clears page, plots file

$$X = \frac{6400i}{k} \quad Y = \frac{6400j}{k} \quad \begin{matrix} \delta j \\ \delta i \end{matrix}$$

C is a cut-off value (Default 4095)

PHOTOMETRY LOAD

n CUTOFF :

FIND (filename)

If windowing required, see [4.1.1]

APPROX

FIT

RESULT

PHOTOMETRY

| | | | | | | |
|-------|-------|------|------|------|------|-------|
| XØ | YØ | PEAK | BASE | HX | HY | R |
| 16.ØØ | 16.ØØ | 3ØØØ | 2ØØ | 8.ØØ | 8.ØØ | Ø.ØØØ |
| Ø.ØI | Ø.Ø3 | I4 | 2 | Ø.23 | Ø.Ø4 | Ø.ØØ7 |

48

XØ, YØ, HX, HY output in pixel units
residuals (related to covariance matrix
 $\sqrt{S^2}$ where S^2 is sum of squared deviations
 $\sqrt{\rho k - 7}$
of previous fit

Loop to [7.2.1] until residuals are sufficiently small.

Creating 'Faked' Gaussians

To create a file containing a digitized gaussian which has been defined by a FIT, immediately after FIT or RESULT.

[7.1]

[7.2]

[7.2.1]

[7.3]

loads relevant programs

Sets cutoff value, default 4095

Selects disk file

Calculates a first approximation

Performs the fit

Outputs results.

4 FITS are usually enough; a single FIT takes about 3 mins. for a 32 x 32 file, 9 mins. for a 64 x 64 file.

If cutoff value to be changed

c CUTOFF !

32767 CUTOFF ! gives c its
maximum value.

If you wish to overwrite the fitted data in a file by its fitted gaussian

FAKE

For every point in the file the
function is calculated and the
value is stored in that point of
the array.

Otherwise create a new file by:-

LISTFILE

Explained in [4.7.2]

n INTER (file 2)

" " "

FIND (file 2) FAKE

Puts gaussian in (file 2)

[7.3.1]

Test gaussians

These can be created by setting up the appropriate parameter values as follows:

Either

FIND (filename)

x₀ RSU X0 !

y₀ UM Y0 !

h_x W0 @ + RSU HX !

h_y W0 l + @ + RSU HY !

)
)
)
) xo, yo, hx, hy in pixels
)
)
)

xo, yo, hx, hy in scaled units described at beginning of [7]

Sets PEAK, BASE, SR

where r = correlation x 1000

Ignore 2nd and 3rd rows of output

Multiplies each value to be stored in the array by $(1 + \frac{\delta(s)}{10000})$ where $\delta(s)$ is normally distributed between $\pm 3s$ with mean zero.

Creates gaussian, stores in disk file.

or
xo XO! yo YØ!

hx HX! hy HY!

p PEAK! b BASE! r SR!

If cutoff to be reset: c CUTOFF!

To check the above entries RESULT

If you wish to add normally distributed noise:

s SD!

Then FAKE

8] Additional Facilities

8.1] Very long scans

For scans with row length too long to fit in the buffer (working storage) the rows must be spilt up into several sections, each section being scanned and stored as a single record on tape (or disk)

Scans as follows:



numbers are record (section) numbers.

Records 0-3, 5-8, 10-13 are of size equal to the buffer size. Records 4, 9 and 14 are shorter since the total row length is in this case, not an integral multiple of the buffer size. (This example might be: 1024 RLENGTH : 3 4500 LWIDE). For scans to tape the scanslot header for records other than record 0 will be incorrect.

```
l BIGSCAN LOAD
  427 LOAD
```

```
K RLENGTH :
```

l (<2,...) is maximum buffer size, loads very-long-scan programs

Buffer size to be used, i.e. standard section length, $K \leq l$. For scans to disk $K \leq 512$.

Sets up parameters as in [3.4] except replace WIDE by LWIDE such that

```
ρ K LWIDE
```

ρ is actual number of rows to be scanned, of K pixels each. LWIDE will automatically divide each row into sections.

For scanning to tape: Scan as in [3.5] except replace WRITE by LWRITE such that

```
LWRITE (identifier)
```

scans onto tape as described above

For scanning to disk: Scan as in [4.7.2] except

```
Replace:      FIND (filename) EQUAL DESTINATION PASS
```

```
by:          FIND (filename) LEQUAL DESTINATION LPASS
```

[8.2] Multiple Scans

Capability to set up for a number of scans to tape as in BIGSCAN, but instead of doing each scan before setting up for the next, store the parameters of each scan in a file called SCANS file, then a single instruction will cause all the scans to be done in turn.

The following parameters may be different for each scan:

The entire "image parameter field": \vec{CN} , \vec{DX} , $\vec{C\phi}$, $\vec{W\phi}$, \vec{WN} , \vec{KIND}

i.e.

(i) size of scan, set by WIDE

(ii) sample intervals, set by DXDY!

(iii) origin, set by either CENTER or c CORNER where c is UR, UL, LL or LR. Scans may be identified by different corners.

(iv) direction of scan, set by VERTICAL or HORIZONTAL

and also (v) The Scansalot 4 ϕ -character identifier

[8.2.1]

K BIGSCAN LOAD
428 LOAD

K is buffer size

CLEAR

Clears scans file

Set up as in [3.4] for first scan

[8.2.2]

n ENTER (identifier)

n is scan no. (1 for first scan)
Enters parameters and identifier in file

Set up PDS for nth scan:

Either

HORIZONTAL

or

VERTICAL

or nothing if previous scan was the same direction

Then

ρ K WIDE

omit if same as previous scan

δx δy DXDY!

" " " " " "

Either

CENTER

After positioning PDS over appropriate point

or

either UR or UL or LL or LR CORNER

" " " " " "

Then Loop to [8.2.2] incrementing n by 1 each time until all scans (≤ 16) have been entered.

If tape is a new tape

REWIND

NTP

SCAN

Then

see [3.5] for explanation

Does each scan in turn, writing onto tape, with the appropriate parameters; with a filemark between each scan and two at the end.

CLEAR

Clears scans file for next user.

[8.3] Array flipping

Capability to flip data in an array on disk

in i : flips each row end for end

in j : flips each column end for end

or in both i and j

[8.3.1]

TAYLOR LOAD

431 OPT LOAD

FIND (file 1) SOURCE FIND (file 2) DESTINATION

Loads flipping option

Allocate file 2 by INTER as in [4.7.2]

Either

ROW-FLIP

or

COL-FLIP

" " " " j

or

RC-FLIP

" " " " i and j

Sets up for flip in i

Then

KEEP PASS

Updates index file and
copies array from (file 1)
to (file 2) flipping on the
way,

To flip an image scanned with origin at one corner to be as if it
was scanned from another, use RECTIFY:

FIND (file 1) SOURCE FIND (file 2) UR DESTINATION RECTIFY

Copies (file 1) to (file 2)
flipping it approximately to
make the scan as if it was
scanned from the UR corner.
Using UL, LL or LR instead of
UR would cause the final file
to have the specified corner as its
origin.

To transform an array scanned in Y to one scanned in X:

FIND (file 1) SOURCE FIND (file 2)
DESTINATION TURN

If (file 1) is A and (file 2) is B, the effect of the command TURN is

A_{ij} becomes B_{ji} where subscripts represent row and column numbers.
Due to the simplicity of the method, the operation is
prohibitively slow for arrays $>(64 \text{ by } 64)$

8.4] Tape file graphic record display

This utility combines SCANSALOT with the display words in SLICING (see
[4.9], [4.12]).

454 LOAD

n SKIPS TWIDE

FIRST

Selects tape file, see [4.12]

Sets scales, plots 1st line.
See [4.9]

The following commands are available, and are described in [4.9]

Reads and plots next line

N

Plots current line

PLOT

AXES

BIG GRID

POINTS

FOLLOW

HISTOGRAM

n m X SIZE

n m Y SIZE

There is one additional command:

SAMPLE

Plots every Nth point of current line
n is set by n SPL! (default 5).
Like PLOT, SAMPLE does not actually
read a record, N must be used for
that first.

[8.5]

Coarse file

Description:

Used most often to compare scans from an electronograph with the corresponding area on a cathode map. Disadvantage is that it is not as accurate as ALIGN, but advantage is that cathode map is only scanned once, and file is reused for different high-resolution scans. ALIGN is generally used to eliminate fine variations within a single scan area while the coarse-file method is used to eliminate grosser variations across a larger area of photocathode.

The whole or a large part of the photocathode map is scanned at a coarse spatial resolution and stored on tape or disk.

A second scan at a higher spatial resolution is stored on disk or tape. The program creates a third file, on disk, containing photocathode map points corresponding to those of second scan, the numbers being the weighted average of the 4 nearest-neighbour points in the coarse file.

8.5.1] Creating the coarse file

Do a normal scan as in [3.4], [3.5] or [4.7.1], [4.7.2] except: ZILCH over a standard mark common to all electronographs to be scanned, (for example, upper left fiducial mark). In addition electronographs must be approximately aligned with one another using the fiducial marks or photocathode defects. The scan must be in X and must be specified by either the CENTER or the UR CORNER.

8.5.2] Creating the nearest-neighbour file from the parameters of a scan (called "source scan") done from a corresponding electronograph.

The source scan is done either onto tape or disk. The sample intervals must be less than those of the coarse file. The source scan must be in X and the origin must be defined by either CENTER or UR CORNER. ZILCH must be done at the same standard point as was done for the coarse file.

The nearest neighbour file may now be created using the data from the coarse file and the parameters of the source scan.

n BIGSCAN LOAD

457 LOAD

If source scan on dish: FIND (filename) SOURCE

Otherwise, i.e. on tape: n SKIPS TWIDE SOURCE

If destination file not allocated in index. use LISTFILE n INTER (file 2)

FIND (file 2) DESTINATION

n is buffer size

loads coarse file routines

Selects source scan file

Selects source scan file see [4.12] for explanation

See [4.7.2] for explanation

Selects destination file

Sets up coarse file parameters

If coarse file on disk: FIND (coarse file) SECOND

Sets up coarse file parameters

If coarse file on tape: n SKIPS TWIDE SECOND

Calculates and stores nearest neighbour file in (file 2)

Then NEAREST

To use (file 2) created by NEAREST to take out effects of photocathode response from (file 1) do TAYLOR LOAD PROCESSING LOAD and see section [4.8]

[8.6]

Density histogram


(Disk files only)

DHIST LOAD


FIND (filename)

DENSITY

Sets up array of 'buckets', goes through file counting number of points in file with densities in range of each bucket, and plots results as a histogram of number against density. The default bucket size is 10 (minimum value is 8, set by n 'QUAN !') i.e. an ordinate of 200 means that there are 200 points between the abscissa density value and that value +10 (or + the value of QUAN)



Example of DENSITY



Example of PLOT

Data may be replotted as follows:-

| | | | |
|-----|-----|-----|------|
| a | b | ROI | PLOT |
| a | b | SEE | PLOT |
| min | max | Y | SIZE |

ENTIRE PLOT

Plots density values a to b only, without change of scale.

Plots density values a to b, expanding scale to fill screen.

Sets Y scale

Equivalent to Ø 4Ø95 SEE PLOT DENSITY plots all values from Ø to 4Ø95

1978 July 26