

S.P.W.


# PSO

USER GUIDE TO PROGRAMS FOR INTERACTIVE IMAGE PROCESSING  
AND MICRODENSITOMETRY ON THE PSO PCS


JONATHAN MCGOWELL

KEITH TAYLOR

JULY 1979



*Handwritten:* 1000  
1000  
1000



## INTRODUCTION

[0.1] Conventions:

- i) User typed FORTH instructions in capitals.
- ii) Terminal response in heavy type.
- iii) Comments and descriptions in lower case.
- iv) In user typed FORTH instructions, lower case symbols such as: n, m, x, y stand for user typed numbers and are defined in accompanying comments. Numbers are to be assumed to be in 16-bit format (integers, unpunctuated) unless stated otherwise.
- v) PDS coordinate directions defined by x and y.
- vi) Size of 2-D data arrays on disk or tape defined by p (for rows) and k (for columns).
- vii) Pixel identification within those arrays given by j, specifying row number and i, specifying pixel in row (column number).
- viii) Some parameters are specified by 2 components. In this case they will be identified using vector notation in the comments. For example  $\overline{MN}$  is a parameter representing the dimensions of a 2-D array or scan. Its components are 2, 16-bit numbers referenced by  $\overline{MN}$  (k - direction)  $\overline{MN}$  l + (p - direction).

[0.2]

Instructions:

FORTH instructions are single words, delimited by spaces. The words will be executed on typing a carriage return (represented here by (cr)). On completion of execution, if execution has not caused any errors, output will be OK. Several instructions may be typed on a line, separated by spaces. They will be executed in turn. FORTH uses a push-down stack for number operations. If an instruction is encountered which is not recognised, the word is echoed back followed by a question mark, and execution is interrupted.

[0.3]

Number entry:

FORTH uses single-precision (16-bit) and double precision (32-bit) signed integers. 16-bit integers are in the range  $-32768 \leq n \leq 32767$ . These are used for most purposes but sometimes, 32-bit numbers are needed, (for example, PDS co-ordinates). 16-bit integers are entered with no punctuation except a leading minus sign where applicable; ie. the same as FORTRAN integers (but there must be no leading plus sign). 32-bit integers are entered by including punctuation either , . - or / . Any number which includes punctuation (except for a single leading

*Jonathan C. M. Powell*

minus sign which means a negative number for both 16- and 32-bit numbers) will be treated as a 32-bit number, but it will be treated as an integer, ie. position of the punctuation will be ignored, hence the numbers 12.34 1.234 12,34 are all taken to be the same number, a double-precision 1234. To avoid confusion, it is probably simplest just to terminate with a point (eg. 1234.)

[Ø.4]

Basic operations:

Arithmetic is reverse polish; integer only.

Typing a number puts it on the stack.

Typing . prints out the top of the stack.

Typing (number) ! stores (number) in (address).

Constants To put value of constant, R on stack just type

```
R
n R !
```

stores n in R

Variables Typing the name of a variable, DX puts address of DX on stack. To put value of DX on stack type

```
DX @
n DX !
```

stores n in DX

The command ? is equivalent to @ . hence it prints out value in address currently on stack; eg. DX ? prints current value of DX:

To set number base, type DECIMAL, OCTAL or HEX.

For more detailed explanation of subjects in [Ø.2] to [Ø.4], see FORTH course notes, FORTH reference manual.

[Ø.5]

FORTH words:

The FORTH compiler only stores the first three letters of each word and its length thus, FILE is treated the same way as FILEL; FORTH cannot tell the difference. FILE1, FILE2, FILE3, FILE4 are treated the same, but FILEAR is not (different length). This is useful when a long word is spelt incorrectly but can lead to confusion. Also note: CENTRE = CENTER DISK = DISC etc. This applies to 'word', ie. FORTH commands, but not to alphanumeric strings. These strings are

usually used as identifiers or file names. They are only valid if used after a command requiring such a string; such commands include FIND, WRITE, CREATE, ENTER. In this documentation, notation of the form name or filename will refer to one of these strings (eg. FIND filename might be FIND GALAXY. These strings should not have any spaces in them, except where mentioned in the documentation.)

[Ø.6] Interrupt:

"Control C" acts as interrupt. This may not work when tape is being used.

[Ø.7] PDS control:

For user to move PDS stage, the SELECT (blue and white) button must be on 'manual'. For the computer to move it, it must be on 'auto' (eg. when doing a scan). Remember to turn control over when using a command which requires stage movement, as the documentation does not mention this. Control can be turned over before or after the relevant command has been sent.

[Ø.8] Images:

The same basic instruction, PASS, is used to transfer scan data from PDS to tape, or from disk to terminal display, etc, ie. from one representation to another. The format of the instruction is (unit1) SOURCE (unit2) DESTINATION PASS to transfer data from (unit1) to (unit2); the format (location, array size, scan sample interval, etc.) is determined by a set of parameters; each unit has its own set. The whole instruction may not always be typed explicitly, eg. in plotting routines the word PLOT may be defined as SOURCE SCREEN-DESTINATION PASS, and so DISK PLOT will transfer the data from the disk to the screen in graphic form.

Each unit name is called an 'image'. A list of 'images' follows.

- Images
- PLATE The PDS stage; PLATE SOURCE (unit2) DESTINATION PASS does a scan
- DISK A file on disk. The word DISK is used implicitly in the instruction FIND filename which also selects a specific file.
- TAPE A file on tape
- TERMINAL The terminal; for digital output
- SCREEN The Tektronix terminal; for graphic output

Also:  
INPUT  
OUTPUT } - Low-level functional images used by PASS  
OPERAND }

Parameters (these can be treated as 'variables')

LOCATION Disk block number for DISK image; not used for other images. Value is  $500 \leq n \leq 2435$ , lower disc or  $5000 \leq n \leq 9500$ , upper scratch disc.

CN Size of scan in pixels; CN contains number of pixels in a row  
CN 1+ contains number of rows in the scan

DX Sample interval in microns:  
DX contains sample interval in X  
DX 1+ contains sample interval in Y

C0 Origin of scan in PDS co-ordinates. This is the first pixel in the array, designated pixel 0, row 0, (for an nxn array pixels and rows are numbered through n-1.) C0 contains 2 double-precision numbers: C0 contains the X co-ordinate. C0 2+ contains the Y co-ordinate.

W0 A 'window' or subset of the array (on disk, tape, etc.) may be defined. The origin (first pixel of first row) in pixels relative to the origin of complete array is contained in W0. W0 contains origin in pixels, W0 1+ contains origin row number. ie. "i", "j" of origin respectively, see [0.1]. W0 and W0 1+ are always positive.

WN Size of window in pixels; WN contains number of pixels/row of window.  
WN 1+ contains number of rows in window.

KIND Specified direction of scan, ie. scan in x or scan in y.

Current Image

The parameters listed above are not ordinary variables because each of the 8 images mentioned above has its own set of parameters. For example, typing DX will reference the DX of only one of these images, called the "current image". An image is made the current image by typing its name (or using a word which itself used an image name). For instance, TAPE DX ? will print out the current DX value for TAPE, while PLATE DX ? will print out that for PLATE, which may be different. Words which change the current image are:

FIND makes DISK the current image  
SEF makes TAPE the current image

PASS makes INPUT the current image

PLOT makes OUTPUT the current image after setting OUTPUT to be equal to SCREEN

SOURCE copies parameters from current image to INPUT and leaves the latter as current image

DESTINATION copies parameters from current image to OUTPUT and leaves the latter as current image

SECOND copies parameters from current image to OPERAND and leaves the latter as current image

USING FORTH

[1] Bootstrap Loading

[1.1] \*\*\*\*\*  
\*\*\*\*\* Follow 'starting up' instructions in black PDS manual up to point where disc  
\*\*\*\*\* is loaded and running: use KTFORTH disc  
\*\*\*\*\*

[1.2] Disk Bootstrap

\*\*\*\*\*  
\*\*\*\*\* Press console switches  
\*\*\*\*\* CTRL HLT/SS (press CTRL and HLT/SS simultaneously)  
\*\*\*\*\* CTRL BOOT  
\*\*\*\*\* Response on VDU screen will be a line of numbers, and then @  
\*\*\*\*\* .....@  
\*\*\*\*\* DK  
\*\*\*\*\* Type  
\*\*\*\*\*

HELLO ?

- ..If
- ..Go to [1.4]
- ..Otherwise go to [1.3]

[1.3] Bootstrap Failure

Hit "Return" key

OK

- ..If
- ..Go to [1.3.1]
- ..Otherwise try reloading the disk;  
see 'closing down' and then 'starting  
up' again in black PDS manual

[1.3.1] Block Buffers Corrupted

33350 4010 ERASE

Should come back OK

Return to [1.2]

[1.4]

Loading

JCM LOAD

.. If this fails, e.g. response is JCM? OCTAL  
go to [1.3]

.. Otherwise output should be

```

* (SYSTEM MESSAGES)
*
System messages
0 (USER MESSAGES)
1
2 User messages

```

09/07/79 22:29 OK

System date & time

.. If you do not wish to correct date  
and time (this is completely optional)  
go to [1.5]

.. Otherwise, for example

```

5/09/78 NOW
18:07 UT

```

Sets date to be 1978 Sep 5  
Sets time to be 18:07; GMT and BST are  
alternatives to UT and are identical,  
computer does not know about time zones

To print out system date, time

```

TODAY @ .DATE 5/09/78
.TIME 28:08

```

\*\*\*\*\*  
User messages are inserted by users to inform other users about possible  
software faults or particular disk files that contain data which other users  
should be especially careful not to overwrite.

To input a line of text to line n of the user message file, (1 ≤ n ≤ 15) use  
the 'Notify Users' command:

```
n NU/ space text / (cr)
```

To input a message to the programmers message file to notify programmers of





Options

\*\*\*\*\*  
Select an option; loading is described under the paragraph number given.  
Subsequently loading another option removes the programs for the first option.\*  
\*\*\*\*\*

<u>Option name</u>	<u>Brief description</u>	<u>Go to</u>
Bigscan	Scan a plate; store data on tape, or on disk Maximum row size is 12,500 pixels Result may be averaged down by factor of $2^n$ in x and y, while scanning.	[ 3 ]
Taylor	Scan a plate, store data on disk Maximum row size is 512 pixels Align a plate, by creating coordinate transform between it and a reference plate. Graphic output; range of additional sub-options are available, eg averaging down in x and y as in Bigscan, adding, subtracting and dividing two images; contour plots, moving arrays between disk and tape, defect removal, disk file graphic display	[ 4 ]
Smoothing	Averaging digitized scan stored on disk, without shrinking array (boxcar convolution)	[ 5 ]
3D	3D plot of a scan stored on disk or tape (several versions)	[ 6 ]
Photometry	Fitting 2D Gaussians to digitized scan (on disk or tape)	[ 7 ]
Additions		[ 8 ]
JPlot	Grey scale plotting routine	[ 9 ]

Error conditions

[ 10 ]

[3] BIGSCAN

[3.1] Loading BIGSCAN

[3.2] Scanning

BIGSCAN LOAD

Wait for the OK before proceeding.

Set up over desired scan

Scan may be specified by center or any corner.

To zero PDS coordinates at this point

ZILCH

To scan in x, go to [3.2.1]

VERTICAL

To scan in y  
Go to [3.2.2]

HORIZONTAL

[3.2.1] If you have most recently scanned

in y, then

Otherwise, "HORIZONTAL" is default  
Continue

[3.2.2] To set PDS velocity;

Either

or

To set scan mode;

Either

or

n MM/S

Sets velocity in mm/s<sup>-1</sup>

m VEL :

Sets velocity in PDS units  
(m = 255 n/60)

ESCAN

Default. Scans in 1 direction only  
increasing x, (or y if VERTICAL)

FSCAN

Scans back-and-forth, and flips rows  
in working storage so final result is  
same

Then set up scan parameters as follows:

PLATE

makes PLATE the current image (see [0.8])

ρ κ WIDE (or ρ κ DIM)

Size of scan; ρ rows with κ pixels in  
each row.

Sets parameters  $\tilde{MN} (=CN)$

Sets sampling interval in microns  
dx is interval in x-coordinate  
dy is interval in y-coordinate  
(sets parameters DX)

Ignore output column marked 'ORIGIN'

6x 6y DXDY

To check the above parameters

STATUS

.. If you wish to define the scan by  
its center

Position over center  
CENTER

.. Then either

.. Or

y. x. \$CENTRE

Uses coordinates y. x. (with the  
decimal point) as scan centre instead  
of actual PDS position. Otherwise  
same as CENTRE

Go to [3.3]

.. Otherwise define the appropriate  
corner

Position over corner

.. Either

UR

Upper right corner; default

.. or

UL

Upper left corner

.. or

LR

Lower right corner

.. or

LL

Lower left corner  
These define corner of scan, as seen  
on PDS viewing screen, not on plate  
itself. They modify the signs of PLATE  
parameters DX

Then either

CORNER

Stores current PDS coordinates in PLATE parameters C0. This position will be the first pixel ( $k = \emptyset$ ;  $p = \emptyset$ ) of the scan.

or

y. x. \$CORNER

AS CORNER but uses y. x. as origin rather than current PDS coordinates

[3.3]

Scanning

For scanning to disk, see the description under TAYLOR LOAD, section [4.7.4]

If tape is a new tape

REWIND

Rewinds tape, initializes variables

NTP

Writes 1 record and 2 filemarks, then backspaces twice. The record is for compatibility with RT tapes which have no filemark before the first file.

Go to [3.3.1]

Otherwise the tape may be positioned anywhere before the final double filemark.

[3.3.1] If you wish to average down during

scan go to [3.4]

Otherwise

WRITE identifier

'identifier' is a string of  $\leq 40$  characters (including spaces). Follow this by carriage return.

(i) Tape will run forward until a double filemark is found, and will then backspace so that the second filemark is overwritten.

(ii) The scan will be performed using the PLATE parameters set previously. The data will be written on to tape in Scansalot format.

(iii) After the scan is completed, a double filemark will be written at the end of the file, and the tape will then be backspaced twice, so that a second WRITE will not require any repositioning of the tape.

BIGSCAN LOAD OK  
ZILGH OK  
3 MM/S OK  
ESCAN OK  
PLATE 64 64 WIDE 10 10 DXDY OK  
CENTER OK  
REWIND OK  
NTP OK  
WRITE TEST SCAN JOM 79RUG4 OK

VERTICAL OK  
5 VEL 1 FSCAN OK  
PLATE 32 28 WIDE 12 5 DXDY OK  
LI 4512 -34128. SCENER OK  
WRITE VERTICAL FSCAN TEST





[4] TAYLOR Options

[4.1] Loading

TAYLOR LOAD

Loads PDS scanning routines, Tektronix graphics routines, least squares fit routines. STARS file routines, and disk file index.  
Max. no. of pixels/scan row is 512.

[4.2] Sub-options

There are a number of additional facilities which may be loaded as well as the programs in TAYLOR LOAD. Only one of these options may be loaded at a time, subsequently loading another one will remove the first.

..Either (i)

PROCESSING LOAD

Loads programs to average in i,j and to add, subtract, divide and normalize two disk files

Go to [4.3]

..or (ii) Either

REMOVAL LOAD

or

SLICING OPT LOAD

See [4.9] before loading either  
Loads programs to display data from a file line by line, and to remove defects by linear interpolation.

Go to [4.9]

..or (iii)

CONTOUR LOAD

Loads programs for contour plot of disk file.

Go to [4.10]  
..or (iv)

SCANSALOT LOAD

Facility to move disk files to tape, and vice versa, read tape files etc.  
Also scans direct to tape.

Go to [4.12]

[4.3]

..If you wish to align a plate before  
.. a scan, so that corresponding arrays  
.. may be obtained from different plates  
.. for the purpose of adding two plates  
.. of the same object, or for the norm-  
.. alization of electronographs using a  
.. photocathode map go to [4.4]

..Otherwise, to scan to disk without  
.. using the plate alignment facility,  
.. go to [4.7]

[4.4]

ALIGN

```
*****  
Description of method: 3 or more points (fiducial marks or stars) are entered  
into a file, called the STARS file, followed by the scan origin coordinates.  
A scan is done. Corresponding points are then entered in file for a second  
plate. Least-squares fitting is used to determine angle of rotation needed  
for the plates to be aligned (the translative component of the coordinate  
transformation is taken care of by the software but the rotation must be made  
zero by physically rotating the stage.) After stage is rotated, re-enter  
points and iterate until angle is acceptably small. (Usually 2 or 3 iterations  
are enough.) Then the corresponding scan is done.
```

```
*****  
(NOTE: the stars file is stored on disk so the data will not be lost until  
they are overwritten by a user.)  
*****
```

Continue

[4.5]

LINK-CLEAR

Cleares out stars file.  
If there is a set of data already in  
it (eg. you are aligning different sets  
of plates) omit this instruction, and  
place the new data later in the file.

Position over the first fiducial mark  
(or star)

If it is a star, you may wish to fit gaussians to it to determine the center point. The word 'mark' will be used (for this section and the next) to mean fiducial mark or star.

.. If gaussians to be fitted

1 ENTER name

'name' is  $\leq 6$  characters long and is for user identification only, not used in programs. This will do a scan 640x1 long (64 pixels) in y (center of scan line is the current PDS position), fit gaussian to the scan line data, plot the data points and the fitted curve on the Tektronix screen; and then

do similar scan in x at y value determined by first fit, fit gaussian, type out co'ords of fitted center and enter them and name into record 1 of STARS file. If you are keeping old data in records 1 through n of stars file, you can start at record n + 1 thus n + 1 ENTER name

.....  
Go to [4.5.1]

.. Otherwise

Center over mark  
1 \$ENTER name

'name' is  $\leq 6$  characters long and is not used in programs. Enters name and current PDS coordinates into record 1 of STARS file. If you are keeping other data in the first n records of STARS file, you can start at record n + 1, thus n + 1 \$ENTER (name)


[4.5.1] To check contents of STARS file at any time

STARS ALL

Proceed to next mark

Follow by two carriage returns. Page will automatically be cleared when return key is pressed for a second time

TAYLOR LOAD OK  
PROCESSING LOAD OK  
LINK-CLEAR OK  
1 ENTER SMALL



Example of ENTER

ROYAL GREENWICH OBSERVATORY

PAGE 1 12/06/78

STARS FILE

#	NAME	X	Y	X	Y
1	FLU		Y		Y
2	FRU	-24778	2	-24767	-11
3	FRL	-24534	-24628	-24607	-24494
4	GALaxy	-9478	-11933	-9547	-11916
5	STAR1	1426	-4209	1356	-4182
6	STAR2	-22551	-5237	-22625	-5210
7	STAR3	-21128	-10525	-21180	-10495
8	H1#1	-10825	-11999	-10991	-11983

These columns not used

OK

Example of STARS ALL





```

: : If gaussian fit required      : OK
: :                               :
: : Go to [4.6.4]                :
: :                               :
: : Otherwise                     : $OK

```

[4.6.4] Then loop back to [4.6.3] until all marks have been stored.

```

Then either
.....
or
Ø i j k ..... 1 NROTATE

```

```

ROTATION
X2      Y2      ANGLE
23.4   -14.2   3.760
12.0   15.0   0.218
137

```

Rotate PDS stage through the appropriate angle

Scans, fits, displays, stores in STARS file (second set of columns). Record number was incremented by NEXT.

Stores current coordinates in STARS file (second set of columns). Record number was incremented by NEXT, so will be correct automatically.

STARS file records n through m inclusive are used to fit a transformation between the two sets of coordinates. If first record was record 1 and 3 marks were entered, 1 3 ROTATE

Fits records i, j, k, ..... 1

Output From ROTATE  
 The number immediately under ANGLE represents the angle or rotation in degrees. X2, Y2 are offsets and second, third rows are residuals. The only number the user needs is the angle.

ave angle = rotate clockwise (plattenn)  
 If angle is large, accuracy in rotation is not important due to the approximation in the transformation.



Iterate by returning to [4.6.1];  
repeat until angle is sufficiently  
small.

Then

NEXT

This will take the stage to the  
appropriate scan position (i.e. the  
position corresponding to the scan  
location entered for the first plate.)  
The second scan may now be done into  
a different disk or tape file. (see  
[4.7] or [4.12]). If several scan  
positions were stored in STARS file for  
the first plate, repeating NEXT will  
move stage to the corresponding position  
so all the scans may be done in turn.

[4.6.5]

Additional ALIGN commands

y. x. GO

Where y, x are 32-bit PDS coordinates;  
moves stage to these coordinates. This  
command is available whenever PDS  
routines are loaded.

n IGO

Moves stage to coordinates given by  
record n of STARS file (1st set of  
coordinates).

n 2GO

Moves stage to coordinates given by  
record n of STARS file (2nd set of  
coordinates).

Graphics mode: The data points for the gaussian fit in ENTER, FIRST, OK are  
plotted out as crosses. ("POINTS MODE") Alternatives are  
listed below. Using the word sets the plotting mode until  
another of the words is used or the graphics are re-loaded.

POINTS

FOLLOW

HISTOGRAM

The commands AXES and BIG GRID will plot axes and a grid on the screen.

[4.7] Scanning

Position PDS over point defining first scan. This point may be the centre or any corner. If ALIGN has been used this will already have been done.

```
*****
* To zero the PDS coordinates at any time use command ZILCH or by
* pressing the appropriate buttons on the console displaying the
* coordinates.
*****
```

To scan in x go to [4.7.1]

To scan in y type

Go to [4.7.2]

[4.7.1] : If you have most recently scanned in y then type

Go to [4.7.2]

Otherwise, scanning in x is default; it is not necessary to type HORIZONTAL unless VERTICAL has been used earlier.

[4.7.2] To set PDS velocity

Either  
or

To set scan mode

Either  
or

VERTICAL

HORIZONTAL

n MM/S

m VEL :

ESCAN

FSCAN

Sets y-scan flag

Clears y-scan flag

Sets velocity in mm s<sup>-1</sup>

Sets velocity in PDS units (m=255n/60)

Default. Scans in l direction only (increasing x, or y if VERTICAL)  
Scans back-and forth, and flips

rows in working storage so final result is same.

Then set up scan parameters as follows

PLATE

Makes PLATE the current image (see [Ø.8]).

$\rho$   $\kappa$  WIDE (or  $\rho$   $\kappa$  YXSIZE)

Size of scan;  $\rho$  rows with  $\kappa$  pixels in each row. Sets parameters  $\overline{WN}$ , and sets  $\overline{CN}$  equal to  $\overline{WN}$ . Note that first parameter is actually no. of rows and not the Y-direction.

$\delta x$   $\delta y$  DXDY

Sets sampling interval in microns.  $\delta x$  is interval in x-coordinate,  $\delta y$  is interval in y-coordinate (sets parameters DX)

To check the above parameters

STATUS

Ignore output column marked 'ORIGIN'

..If you wish to define the scan by its center  
..Then either

Position over center  
CENTER

Stores origin (ie. coordinates of most -ve corner) by setting parameters  $\overline{CØ}$ , using current PDS position (x, y) and  $\overline{WN}$  and  $\overline{DX}$ ,  
Such that  $x_0 = x - \frac{1}{2}w_x \delta x$   
 $y_0 = y - \frac{1}{2}w_y \delta y$   
where ( $w_x, w_y$ ) are ( $\kappa, \rho$ ) or ( $\rho, \kappa$ ) depending on scan direction.

..or

y. x. \$CENTRE

As CENTRE but uses coordinates y. x. instead of current PDS position

Go to [4.7.3]

..Otherwise define the appropriate corner

Position over corner

..Either  
...or  
...or  
...or

UR  
UL  
LR  
LL

Upper right corner: default  
Upper left corner  
Lower right corner  
Lower left corner

ROYAL GREENWICH OBSERVATORY

PAGE 1 12/06/78

IMAGES INDEX

# NAME	LOC'N	ORIGIN	X	Y	DIR.	CORN.	SIZE	X	Y	+X	+Y
1 CHRIS	435	0	0	0	H	UR	128	128	10	10	10
2 DRETT	516	648	648	648	H	UR	64	64	60	60	60
3 COPY	520	-1920	-1920	-1920	H	UR	64	64	60	60	60
4 JCM	527	0	0	0	H	UR	128	128	48	48	48
5 LONE	470	0	0	0	U	UR	512	4	10	10	10
6 MED21	477	-320	-320	-320	H	UL	32	32	10	10	10
7 MISS#1	458	-1160	1160	1160	H	UR	59	50	99	99	99
8 PK	432	-486	-521	-521	H	UR	64	64	10	10	10
9 TEST	461	14	-1	-1	H	LR	64	64	10	10	10
10 TWIN	458	328	-328	-328	H	LR	64	64	10	10	10
	432	-1000	1000	1000	H	UL	32	32	10	10	10
	432	-1160	1160	1160	H	UL	32	32	10	10	10
	439	0	0	0	U	UR	64	64	10	10	10
	439	320	320	320	U	UR	64	64	10	10	10

OK

(Name)

Starting block  
ending block

X origin  
X center  
Y origin  
Y center

Horizontal or  
Vertical scan?

Which corner?  
(corner=UR)

X pixels  
(K)

Y pixels  
(P)

6x  
microns

Example of LISTICE



To remove old files in the index  
that are no longer used

REMOVE filename

n is starting block number (500 ≤ n ≤ 2435). Block nos. 5000 to 9500 access the upper (scratch) disc. filename is a string of ≤ 8 characters. This creates a file accessed by the command FIND filename (Do not begin filename with an S!)

Remove entry from index; Files should be removed as soon as they are finished with. This command can be used at any time.

n m LISTB

Restricted LISTFILE that only lists files occupying disk space between blocks n and m. (n < m)

n UP LISTB

Lists files occupying space from block n upwards.

..If you have done a PROCESSING LOAD  
..(or do one at this point) and wish  
..to use the averaging capability  
..for this scan, go to [4.7.5]

..Otherwise

FIND filename EQUAL DESTINATION

FIND filename selects the file (the one just created, usually) EQUAL copies INPUT's parameters to DISK and DESTINATION specifies DISK as the OUTPUT image, copying parameters from DISK to OUTPUT.

PASS

Passes INPUT to OUTPUT. Since INPUT is PLATE and OUTPUT is DISK, PASS will do a scan according to the parameters set up, and store the data on disk. Remember to select AUTO on PDS.

..If you are using ALIGN and now  
..wish to align second plate, go  
..to [4.6]

..Otherwise, continue





..If you are using ALIGN and now  
..wish to align second plate, go  
..to [4.6]

..Otherwise, continue

..If you wish to do further scans, go  
..to [4.7]

..Otherwise return to [2] or finish

\*\*\*\*\*  
\* For further capabilities in PROCESSING LOAD see [4.8] \*  
\*\*\*\*\*

[4.8] Further capabilities in PROCESSING LOAD

[4.8.1] Combining disk images

\*\*\*\*\*  
\* For two disk files (File 1) and (File 2) to be arithmetically combined to form  
\* a third file (File 3) use the following instructions. \*  
\*\*\*\*\*

FIND file 1 SOURCE file 1 is identified with INPUT

FIND file 3 DESTINATION file 3 is identified with OUTPUT

FIND file 2 ADD (SUBTRACT) (DIVIDE)  
file 2 is identified with OPERAND by  
using one of the words.

PASS  
PASS combines data in INPUT with data  
in OPERAND and places the result in  
file 3

\*\*\*\*\*  
\* The commands defining the OPERAND are: \*  
\*\*\*\*\*

ADD giving (File 3) = (File 1) + (File 2)  
SUBTRACT giving (File 3) = (File 1) - (File 2)  
DIVIDE giving (File 3) = NUM x (File 1) / (File 3)

```

*****
where NUM is a variable (default 1000) which may be set by n NUM !
*****
Each element in the INPUT array is combined with the corresponding element
in OPERAND array row by row in working storage. The INPUT and OPERAND files
are unaltered.
*****

```

```

.....
..If you wish to overwrite file 1
..with the result then
..
..
..
..
..Otherwise, finish. However, if you
..wish to overwrite (file 2) with the
..result, then

```

```

      FIND file 1 SOURCE'DESTINATION
      FIND file 2 ADD (SUBTRACT, DIVIDE)
      PASS
      FIND file 1 SOURCE
      FIND file 2 DESTINATION
      ADD (SUBTRACT, DIVIDE) PASS

```

[4.8.2] Sum and average of an array

```

*****
SBAR will calculate the average value of an array and both print it out and
store it in the variable NUM. The sum of the array is stored in the double
precision (32-bit) variable XSUM, to see value; XSUM 2@ D.
To subtract a constant background value n throughout the calculation (while
not altering the array itself) use the command
      n BKG !
before using SBAR; default BKG is zero, after changing it, it will stay at
current value until altered again or until PROCESSING is reloaded.
*****

```

```

      FIND filename SBAR AVERAGE VALUE
      = 1025
      calculates sum and average

```

[4.8.2.1] Subtraction of background

To subtract a constant from an array on disk

FINDD filename  
n SUBBK

Subtracts n from array 'filename'

...  
To subtract the average value, (provided the value of BKG is zero)

FINDD filename SBAR NUM @ SUBBK

[4.8.3] Normalization of arrays

\*\*\*\*\*  
For example, for removing variations due to photocathode sensitivity on an electronograph, divide scan array by scan photocathode map and multiply through by average of latter scan.

\*\*\*\*\*  
Suppose first scan is in a disk file called STAR, the photocathode map in a file called SKY, and you wish to place the result in a file called FINAL, which you first need to create with size equal to those of STAR and SKY, starting at block n.  
\*\*\*\*\*

LISTFILE

See [4.7.4] for details.

FINDD STAR

Makes disk current image with STAR's parameters.

n CREATE FINAL

Enters FINAL in index file with parameters equal to STAR.

FINDD SKY SBAR AVERAGE VALUE = 952

Stores average value of SKY by NUM.

FINDD STAR SOURCE FINDD FINAL  
DESTINATION  
FINDD SKY DIVIDE PASS

Divides STAR by SKY.  
Multiplies through by NUM and puts the result in FINAL.

To do further DIVIDE's with NUM at the default value

1000 NUM :





Position cursor on R.H.S. of defect  
Tap key

Program does linear interpolation  
between 2 points.

New interpolated line is stored in the  
disk file and replotted on screen.  
Old line is lost.

\*\*\*\*\*  
\*\*\*\*\* These are all the available words in SLICING \*\*\*\*\*  
\*\*\*\*\* If you have loaded REMOVAL continue to [4.9.2] \*\*\*\*\*  
\*\*\*\*\*

#### [4.9.2] REMOVAL

\*\*\*\*\*  
\*\*\*\*\* This assumes that source of disk file is still on PDS and coordinate system \*\*\*\*\*  
\*\*\*\*\* has not been changed (except for section [4.9.2.4]) \*\*\*\*\*  
\*\*\*\*\*

#### [4.9.2.1]

..... If you wish to enter removed-defects  
..... co-ordinate into DEFECTS file to  
..... initialize that file then

Plot lines until possible defect  
appears

If desired clear page and PLOT

?DEFECTS (cr) (cr)

Displays cursor and contents of  
DEFECTS

LINK-CLEAR

Initializes file

..... Continue

NOX

Puts up cursor

..... Otherwise

Centre cursor over apparent defect

This position is not used in any  
subsequent defect removal and hence  
position of cursor is not critical.

Tap any key (except RETURN)

PDS will move to co-ordinates  
corresponding to that point and the  
user may decide from visual appearance  
of area whether removal of feature is  
desired.

.. If defect to be removed go to  
[4.9.2.2]

.. Otherwise go back to [4.9.2.1] and  
locate another defect candidate

[4.9.2.2] If defect to be entered into  
DEFECTS file go to [4.9.2.3]

.. Otherwise use FIX as described  
before in [4.9.1.1]

[4.9.2.3]

NIX

Identical to FIX except, stores  
current PDS co-ordinates in DEFECTS  
file. (Do not move PDS between NOX  
and NIX otherwise wrong co-ordinates  
will be entered.)

To view DEFECTS file

?DEFECTS

Defects are entered sequentially  
numbered after last defect in file.

To remove further defects go to  
[4.9.2.1]

.. If a second plate is then set up  
.. using ALIGN so that the co-ordinate  
.. translation has been calculated  
.. using ROTATE

n @DEFECT

This will move PDS to position  
corresponding to coordinates in  
record n of defects file, with the  
translation added.  
In this way the user can see if the  
same defects are present on different  
plates.  
If no ROTATE has been performed since  
loading TAYLOR then this instruction  
will take PDS to nth defect of existing  
plate.

.. Otherwise, continue

[4.9.2.4] Fitting a Gaussian to a Scan Line

Plot line by using

N (or PLOT)

APPROX FIT

Calculates 1st approximation and does least squares fit.

CURVE

Plots fitted curve.

RESULT

CENTERING

CENTER	PEAK	BASE	HW
1.52	3000	257	5.42
0.00	0	0	0.00
0			

Outputs values of parameters X0 (value in column CENTER) is distance from center of curve in units of 1/10  $\mu$  (ie. if result is x, pixel number in line is  $x/100x + k/2$  where k is pixel length of row)

```

*****
PEAK  Is height of the peak above the base.
BASE  Is the base. Units for these are the same as the scan data
      (0  $\leq$  density  $\leq$  4095).
HW    Halfwidth, units 1/10  $\mu$ . Second line gives residuals
      Repeat FIT, CURVE, RESULT as required, for example:
      FIT CURVE FIT FIT CURVE FIT RESULT CURVE etc.
*****
  
```

[4.10]

CONTOUR PLOTS

CONTOUR LOAD

FIND filename

Selects disk image to be plotted.

- .. If only a part (window) of the file
- .. is to be plotted, see [4.11] to do
- .. this
- .. Otherwise, continue

.. Either

XWIDTH

Sets size of picture

.. or

SCREEN n WIDTH

Where n is pixel width of screen

n SPACE :

Optional: default 100  
Sets contour interval to be n (ie. difference of n in the numbers in the



array.)

Optional: default - 32768 32767 CLIP  
"low" and "high" are numbers defining  
cutoff points, for example 200 4000  
CLIP will treat all numbers < 200 as  
200, and all > 4000 as 4000.

Plots contour map

low high CLIP

DISK PLOT

[4.10.2] Contour plotting directly from tape

512 MAXIMUM SCANSALOT LOAD  
CONTOUR LOAD

~~n SKIPS TAPES~~  
**SET**  
n SKIPS ~~TAPES~~

This unloads TAYLOR so if you  
subsequently wish to use programs in  
TAYLOR you must repeat TAYLOR LOAD.  
Select tape file; described in [4.12]

..If windowing is desired, see  
..[4.11] , then go to [4.10.3]  
..Otherwise continue

[4.10.3]

SCREEN n WIDTH

n SPACE ! low high CLIP


TAPE PLOT

Described in [4.10.1]

" " "


Plots contour map.

Q 1200 CLIP 100 SPACE 1 DISK PLOT



OK

Example of a PLOT using CONTOUR 42



Given a 2D array  $A_{ij}$  of pixel dimension  $\rho \times \kappa$  (such that  $0 \leq i < \kappa$  and  $0 \leq j < \rho$ ) a subset  $B_{ij}$  of pixel dimension  $a \times b$  (such that  $m \leq i < (m+a)$  and  $n \leq j < (n+b)$ ) may be accessed, provided  $m$  and  $n \geq 0$ ,  $(m+a) < \kappa$  and  $(n+b) < \rho$ .

This is used for plotting part of a large array, combining part of an array with a smaller array fitting a function to part of an array, etc.

Windowing words are:

- n m ORIG
- b a WM
- FULL

Sets  $\tilde{W}$ : defines window origin at row  $n$ , pixel  $m$ .

Sets  $\tilde{M}$ ; defines size of window to be  $b$  rows of  $a$  pixels each.

Sets  $\tilde{W}$  to  $(0,0)$  and  $\tilde{M} = \tilde{CN}$  ie: makes the window the default value, the whole array.

58	53	60	61	62	63	64	65	66	67	68
531	541	554	554	555	554	551	541	551	564	540
532	554	557	552	562	559	570	575	563	565	559
533	567	576	573	591	592	586	587	588	580	584
534	586	592	603	647	635	611	609	607	602	592
535	595	620	642	655	635	655	654	647	632	620
536	624	630	634	691	710	717	717	695	676	655
537	672	700	725	742	753	775	776	764	738	701
538	672	725	789	823	845	836	839	830	794	744
539	721	757	789	823	845	836	839	830	794	744
540	759	810	858	917	946	945	933	908	860	797
541	803	884	956	1000	1077	1100	1047	995	937	858
542	869	936	1077	1177	1245	1305	1364	1194	937	858
543	945	976	1135	1310	1430	1440	1353	1320	1015	960
544	993	1062	1155	1470	1591	1581	1473	1308	1133	988
545	1042	1116	1293	1532	1653	1639	1516	1326	1154	996
546	1090	1152	1345	1572	1693	1614	1473	1295	1118	980
547	1138	1199	1393	1600	1649	1614	1456	1211	1054	939
548	1186	1252	1442	1637	1406	1456	1376	1093	1054	877
549	1234	1307	1490	1675	1442	1276	1203	1093	979	820
550	1282	1362	1538	1713	1485	1093	1053	996	914	760
551	1330	1417	1586	1751	1523	992	927	894	831	700
552	1378	1472	1634	1789	1561	992	850	795	759	710
553	1426	1527	1682	1827	1600	992	779	729	702	671
554	1474	1582	1730	1865	1638	992	730	681	666	644
555	1522	1637	1778	1903	1676	992	682	644	609	612
556	1570	1692	1826	1941	1715	992	641	611	609	582
557	1618	1747	1874	1979	1753	992	600	579	570	549
558	1666	1802	1922	2017	1791	992	559	565	557	547
559	1714	1857	1970	2055	1830	992	518	549	546	526
560	1762	1912	2018	2093	1868	992	477	527	532	526
561	1810	1967	2066	2131	1906	992	436	511	517	519
562	1858	2022	2114	2169	1945	992	395	511	509	505
563	1906	2077	2162	2207	1983	992	354	511	509	505
564	1954	2132	2210	2245	2021	992	313	511	509	505
565	2002	2187	2258	2283	2060	992	272	511	509	505
566	2050	2242	2306	2321	2098	992	231	511	509	505
567	2098	2297	2354	2359	2136	992	190	511	509	505
568	2146	2352	2402	2397	2175	992	149	511	509	505
569	2194	2407	2450	2435	2213	992	108	511	509	505
570	2242	2462	2498	2473	2251	992	67	511	509	505
571	2290	2517	2546	2511	2290	992	26	511	509	505
572	2338	2572	2594	2549	2328	992	15	511	509	505
573	2386	2627	2642	2587	2366	992	14	511	509	505
574	2434	2682	2690	2625	2405	992	13	511	509	505
575	2482	2737	2738	2663	2443	992	12	511	509	505
576	2530	2792	2786	2701	2482	992	11	511	509	505
577	2578	2847	2834	2739	2520	992	10	511	509	505
578	2626	2902	2882	2777	2559	992	9	511	509	505
579	2674	2957	2930	2815	2597	992	8	511	509	505
580	2722	3012	2978	2853	2636	992	7	511	509	505
581	2770	3067	3026	2891	2674	992	6	511	509	505
582	2818	3122	3074	2929	2713	992	5	511	509	505
583	2866	3177	3122	2967	2751	992	4	511	509	505
584	2914	3232	3170	3005	2790	992	3	511	509	505
585	2962	3287	3218	3043	2828	992	2	511	509	505
586	3010	3342	3266	3081	2866	992	1	511	509	505
587	3058	3397	3314	3119	2905	992	0	511	509	505
588	3106	3452	3362	3157	2943	992	0	511	509	505
589	3154	3507	3410	3195	2982	992	0	511	509	505
590	3202	3562	3458	3233	3020	992	0	511	509	505
591	3250	3617	3506	3271	3059	992	0	511	509	505
592	3298	3672	3554	3309	3097	992	0	511	509	505
593	3346	3727	3602	3347	3136	992	0	511	509	505
594	3394	3782	3650	3385	3174	992	0	511	509	505
595	3442	3837	3698	3423	3213	992	0	511	509	505
596	3490	3892	3746	3461	3251	992	0	511	509	505
597	3538	3947	3794	3499	3290	992	0	511	509	505
598	3586	4002	3842	3537	3328	992	0	511	509	505
599	3634	4057	3890	3575	3366	992	0	511	509	505
600	3682	4112	3938	3613	3405	992	0	511	509	505

OK

Example of VCLLA

[4.11.2] Digital Display

FNND filename n n VOILA

Clears page and displays window of file with origin at pixel n, row m. Resets ~~W~~, ~~M~~ to previous values at end

[4.11.3] Examples of Copying Files

FNND (file 1) SOURCE  
FNND (file 2) DESTINATION PASS

Copies (file 1) to (file 2)

FNND (file 1) 32 32 W W SOURCE  
FNND (file 2) 32 32 W W DESTINATION  
PASS

Copies "top left" (1st 32 pixels of 1st 32 rows of file 1) to corresponding area of file 2.

[4.12] SCANSALOT LOAD

\*\*\*\*\*  
\* These commands are also available in BIGSCAN (See [3]) \*  
\*\*\*\*\*

[4.12.1] To set up parameters for a scan to tape, go to [4.7.1]

.....If a tape is a new tape (no files on it to be kept)

REWIND

Rewinds tape, and initialize tape variables.

NTP

Writes 1 dummy record and 2 filemarks then backspaces twice. The record is for compatibility with RTSCAN tapes which have no filemark before the first file.

.....Go to [4.12.2]

.....Otherwise ensure tape is positioned at some point before the final double filemark. Program will look for a double filemark before starting to write.

[4.12.2] Then

WRITE identifier (cr)

Identifier is a string of ≤ 40 characters (including spaces). Tape will run forward until a double filemark is found, and will then backspace so that the second filemark will be overwritten. The scan will be performed using the PLATE parameters set previously. The data will be written onto tape in Scansalot format.

After the scan is completed, a double filemark will be written at the end of the file, and the tape will then be backspaced twice, so that a second WRITE will not require any repositioning of the tape.

\*\*\*\*\*  
\* The rest of this section [4.12] can be used without doing a TAYLOR LOAD, by  
\* doing a 512 MAXIMUM SCANSALOT R-X instead.  
\*\*\*\*\*

[4.12.3] Disk-to tape pass

FNND filename

Selects appropriate disk file

WRITE identifier

Looks for end of tape (double filemark) backspaces once and copies disk file onto tape. Finishes by writing a new double filemark and backspaces twice.

Identifier is ≤ 40 character string including spaces terminated by carriage return.

[4.12.4] Tape positioning (word definitions)

n POS

Immediately after loading, the computer does not know the current position of the tape, so if the tape is not at the beginning, you must tell it which file number it is at, using the command POS.

T-INDEX  
 0 DUMMY FILE ; BEGINNING OF TAPE  
 1 JCM TEST SCAN 1978JUL26  
 2 VERTICAL SCAN, ODD=1  
 3 ANOTHER VERTICAL SCAN  
 4 SIMPLE SCAN  
 5 REWRITE TEST  
 END OF FILES OK

2	1	0	0	0	0	0
64	64	10	10	10-26860.	-3398.	0.
20	32	-12	-5	-14010.	-11834.	0.
20	32	-12	-5	-14010.	-11834.	0.
64	64	10	10	10-26860.	-3398.	0.
64	64	10	10	10	0.	0.

CUN CUN4T 0X 0X4T C3 C0 2T

Example of T-INDEX

REWIND

Rewinds tape, initializes tape variables EOFs and REC.

n SKIPS

n > 0 passes forward over n filemarks. Remains positioned immediately after the last (ie. at the beginning of a file).

n < 0 goes backward over n filemarks then forward 1 space, so positioning tape at the beginning of a file.

-1 SKIPS returns to beginning of current file. (Sometimes!)

SET

To be used after SKIPS (or FILE).

Returns to beginning of current file, goes through file checking number of records and coordinates in header, uses these to set TAPE parameters

WN, CN, C0, DX. Types first identifier returns to beginning of the file.

n SKIPS SET is the taping equivalent of FIND filename for disk files.

T-INDEX

Prints directory of files on tape; rewinds tape and goes forward until double filemark found, then backspaces twice.

n FILE

Positions at end of first record of file n, n is absolute file number with dummy header file being zero. However in some circumstances when interrogating individual records in a file, one goes past the end of file. Subsequent use of FILE may give the wrong result.

LABEL

Reads and prints next records identifier

n RECS

Moves forward n-1 records and reads next record.



[4.12.5]

Tape to disk pass

SHOW

n SKIPS LABEL

- 1 RECS reads next record
- 2 RECS reads next record but 1
- 0 RECS reads record just read
- 1 RECS reads record before one just read.

Types identifier and dumps scan data currently in working storage. To be used after RECS. (N.B. WN may be incorrect if a TWIDE has not been done on file).

Moves n filemarks then prints identifier of first record (to check if file is correct one).


```

*****
***** If disk array not yet allocated, use LISTFILE or n m LISTB to examine images
***** index and then create file by n CREATE filename, as explained in [4.7.2]
*****

```

TDPASS filename

Copies tape file to disk file 'filename'.



[5] SMOOTHING (boxcar convolution)

[5.1]

Averaging in  $i$  and  $j$  by constant factor  $a$ ; disk arrays only.

Description:

Let contents of pixel  $i$ , row  $j$  in array  $A$  be  $A_{ij}$

$A$  is source array,  $B$  is smoothed array

Each array consists of  $\rho$  rows of  $k$  pixels each

$[X]$  means integer part of  $X$

Action of program:

$B_{ij} = A_{ij}$  for  $(i) \quad \emptyset \leq j < [a/2], \emptyset \leq i < k$

$(ii) \quad \emptyset \leq j < [a/2], \emptyset \leq j < \rho$

$(iii) \quad k - [(a+1)/2] + 1 \leq i < k, \emptyset \leq j < \rho$

Partially smoothed s.t.  $B_{ij} = (a(\sum_{\Sigma}^{p-1} (i + [(a+1)/2] - 1) - 1) \sum_{\Sigma}^{k-1} A_{k\ell})^{-1}$   
 $\ell = j - [a/2]; k = [i - a/2]$

for  $[a/2] \leq i < k - [(a+1)/2] + 1; \rho - [(a+1)/2] + 1 \leq j < \rho$

Smoothed s.t.  $B_{ij} = a^{-2} \sum_{\Sigma}^{j + [(a+1)/2] - 1} \sum_{\Sigma}^{i + [(a+1)/2] - 1} A_{k\ell}$   
 $\ell = j - [a/2]; k = i - [a/2]$

for  $[a/2] \leq i < k - [(a+1)/2] + 1, [a/2] \leq j < \rho - [(a+1)/2]$

Inside fully smoothed region:

For odd values of  $a$  -  $B_{ij}$  is the average of the  $a \times a$  square of points in array  $A$  with center  $(i, j)$

For even values of  $a$  -  $B_{ij}$  is the average of the  $a \times a$  square of points in array  $A$  with center  $(i - \frac{1}{2}, j - \frac{1}{2})$  so for even  $a$

the image is translated by  $\frac{1}{2}$  pixel in  $i$  and  $j$ . If an array is to be repeatedly smoothed by an even factor it is recommended that the array be 'flipped' each time before smoothing so that:  $A_{ij} \rightarrow A_{(k-i-1, p-j-1)}$ . Thus the successive translations will cancel instead of adding. How to achieve this is described below.

[5.2] Tape Files

```
*****
This program only works on disk files. To smooth a tape file first copy it on
to disk, either 512 MAXIMUM SCANSALOT LOAD or TAYLOR LOAD SCANSALOT LOAD, and
then go to [4.12.5] for instructions on tape-to-disk pass.
*****
```

[5.3] Loading

```
a SMOOTHING LOAD
a is the full width of the boxcar in
p and k ( $a \geq 2$ ).
```

[5.4] Disk Files

Allocate all disk files you will require using

```
LISTFILE
FIND file 1
n CREATE filename
```

Lists images index; for explanation see [4.7.2]  
 Sets current image parameters to those of file 1. file 1 is the source file ie. the one you wish to smooth.  
 n is block number, select as described in [4.7.2]. CREATE allocates a file called filename the same size as file 1. Use CREATE to allocate all the files you will need.

[5.5.1]

Smoothing

(i) Single smoothing

FINDD file 1 SOURCE

FINDD file 2 CONVOLUTE

FINDD file 1 \$CONVOLVE

(ii) Smoothing into same file

FINDD file 1 \$CONVOLVE \$CONVOLVE  
\$CONVOLVE

(iii)

Flipping and then smoothing, for use with repeated smoothing by even values.

FINDD file 1 SOURCE

FINDD file 2 RCONVOLUTE

SOURCE FINDD file 1 RCONVOLUTE

SOURCE FINDD file 2 RCONVOLUTE

SOURCE FINDD file 1 RCONVOLUTE

SOURCE FINDD file 2 RCONVOLUTE

Defines file 1 as the INPUT file, ie. file to be smoothed.

Copies file 1 to file 2 and smooths file 2. CONVOLUTE is equivalent to DESTINATION PASS \$CONVOLVE (see below).

Smooths file 1, overwriting original data. If you wish to save original data, copy it into another file or use CONVOLUTE as in (i) and then FINDD file \$CONVOLVE.

Repeatedly smooths file 1

file 1 will be overwritten later, to save original first copy.it.

Copies file 1 to file 2 reflecting it on the way so that pixel (i,j) of file 1 is put in pixel (k-i-1, p-j-i) of file 2 ie. reflected in the center point of the array.  
Then smooths file 2

Flips file 2 into file 1 and smooths file 1

Repeatedly smooths array, with the 1/2-pixel translation caused by each smoothing cancelling.





To set the X scale

n m X SIZE

n is the abscissa of the LHS of the screen.  
m is the abscissa of the RHS.

Thus

Ø WN @ X SIZE

will cause a plotted row to fill the screen from left to right (WN @ gives value of k).

n m Y SIZE

set the Y scale.

n m Z SIZE

controls spacing between rows.

Default values:

Ø 32 X SIZE

Ø 8192 Y SIZE

- 64 64 Z SIZE

The pitch and yaw of the image can also be altered.

Pitch inclines the Z axis from the horizontal.  
Yaw inclines the Y axis from the vertical, in the plane normal to Z.  
Default pitch is 60°, yaw is zero.

[ 6.2 ]

Commands

Commands for the three variants are identical except in the case of PITCH and YAW.

SCREEN p k WIDE

No more than p rows and no more than k pixels/row will be plotted.

FIND filename


Selects disk file to be plotted.

..Either for a disk file  
..or, for a tape file


n SKIPS SET

n is relative file number;  
see [4.12] for explanation.





DK



2 PUT usj 30/5 58

[6.2.1] Then for Disc or Tape

..If only part of file to be plotted

n m ORIG

b a MW

see [4.11] for explanation.  
see [4.11] for explanation.

Continue

n m X SIZE

Sets up X axis co-ordinates.

n = abscissa of L.H.S.

m = abscissa of R.H.S.

(Default 0 32 X SIZE)

..Otherwise, continue

n m Y SIZE

(Default 0 8192 Y SIZE)

n m Z SIZE

(Default -64 64 Z SIZE)  
For explanation see [6.1]

..If 3D variant being used

s c PITCH

Default 866 500 PITCH

s c YAW

Default 0 1000 YAW

s = 1000 sin  $\theta$  c = 1000 cos  $\theta$

Go to [6.2.2]

..Otherwise (ie. 3DH or 3DJ variants).

nn.nn PITCH

Default 60.00 PITCH,  
nn.nn is angle in degrees. There must be exactly 2 decimal places - do NOT omit trailing zeros.

nn.nn YAW

Default 00.00 YAW

[6.2.2] ..If a disk file then

PLOT

Clears page, plots file.

Go to [6.2.3]

..Otherwise for a tape file

PLOT

Clears page, plots file.

[6.2.3]

To replot using different window, scales, pitch or yaw: go to [6.2.1]

Description: Least squares fitting of a 2-D gaussian to scan data from a disk file. The file may be windowed (see [4.11]).

Consider a 2-D array of data  $A_{ij}$  where  $0 \leq i < \kappa$  and  $0 \leq j < \rho$

Fitting Parameters

$X\emptyset = 6400i/\kappa$

$Y\emptyset = 6400j/\kappa \cdot \delta j/\delta i$

where ( $X\emptyset$ ,  $Y\emptyset$ ) is centre of gaussian and  $\delta j/\delta i =$  ratio of sample intervals in microns of the  $j$  and  $i$  direction.

PEAK = Height of peak above base (units as in  $A_{ij}$ )

BASE = Base (units as in  $A_{ij}$ )

HX = Half width at half height in  $i$  direction (units as in  $X\emptyset$ )

HY = Half width at half height in  $j$  direction (units as in  $Y\emptyset$ )

SR = Correlation. Stored as integer. Value multiplied by 1000.

Function:

$$f(i, j) = (\text{BASE}) + (\text{PEAK}) \cdot 2^{\xi} \leq C$$

$$\text{where; } \xi = \left[ \frac{X - X\emptyset}{HX} \right]^2 + \frac{SR(X - X\emptyset)(Y - Y\emptyset)}{HX \cdot HY} + \left[ \frac{Y - Y\emptyset}{HY} \right]^2$$

$$X = 6400i/\kappa ; \quad Y = 6400j/\kappa \cdot \delta j/\delta i$$

C = cut-off value. (Default 4095)

[7.1] To initiate PHOTOMETRY

PHOTOMETRY LOAD

Loads relevant programs

[7.2] n CUTOFF !

n CUTOFF !

Sets n to be cutoff value. Default:  
n = 4095

..

FIND filename

Selects disk file

..If windowing required; see [4.11],  
: then go to [7.2.1]  
: :  
: Otherwise, continue

[7.2.1]

APPROX

Calculates a first approximation

[7.2.2]

FIT

Performs the fit

RESULT	PHOTOMETRY
X0	Y0 PEAK BASE HX HY SR
16.00	16.00 3000 200 8.00 8.00 0.000
0.01	0.03 14 2 0.25 0.04 0.007
48	

Outputs the result

X0, Y0, HX, HY output in pixel units.  
Residuals; related to covariance matrix  
VS2 where S2 is sum of squared deviation  
/ok-7 of previous fit.

Loop to [7.2.2] until residuals  
are sufficiently small.

4 FITS are usually enough; a single  
FIT takes about 3 mins. for a 32 x 32  
file, 9 mins. for a 64 x 64 file.

[7.3]

Creating "Faked" Gaussians

```

*****
* To create a file containing a digitized gaussian which has been defined by a
* FIT, immediately after FIT or RESULT
*****

```

To change cutoff value

c CUTOFF !

32767 CUTOFF ! gives c its maximum  
value.



To add normally distributed noise

s SD :

Multiples each value to be stored in the array by [1 +  $\delta(s)/1000$ ] where  $\delta(s)$  is normally distributed between  $\pm 3s$  with mean zero.

Then

FAKE

Creates gaussian, stores in disk file.

[7.4]

Graphic Display of Fit

```

*****
After fit is complete, you may use all the commands in SLICING (except FIX)
to display the file, and also the command CURVE which plots the fitted line
corresponding to the current line plotted. The command XC goes through the
whole file doing FIRST CURVE N CURVE N ..... etc, clearing the page between
each row, thus allowing you to view the goodness of fit of the whole file.
*****

```

[8] Additional Facilities

[8.1] Multiple scans

Capability to set up for a number of scans to tape as in BIGSCAN, but instead of doing each scan before setting up for the next, store the parameters of each scan in a file called SCANS file. A single instruction will then cause all the scans to be done in turn.

Max. no of scans allowed is 79.  
The following parameters may be different for each scan:

The entire "image parameter field": CN, DX, CO, WO, WT, KIND

ie. (i) size of scan, set by WIDE (or DIM)

(ii) sample intervals, set by DXDY !

(iii) origin, set by either CENTER or c CORNER where c is UR, UL, LL or LR. Scans may be identified by different corners.

(iv) direction of scan, set by VERTICAL or HORIZONTAL.

and also

(v) the Scanslot 40-character identifier

(iv) whether the scan is averaged or not, and by how much

[8.1.1]

BIGSCAN LOAD 428 LOAD

CLEAR

clears scans file

Go to [3.2] to  
Set up PDS for first scan  
(ie. MM/S, WIDE, etc.)



[8.1.2]..Either

n ENTER identifier

n is scan number (1 for first scan)  
Enters parameters and identifier in  
file.

.....  
or

a n A-ENTER

a is value to be averaged by  
(= 2<sup>n</sup>) in x and y

Set up PDS for (n + 1)<sup>th</sup> scan

Then

.....  
or

HORIZONTAL

Default

.....  
or

VERTICAL

.....  
or nothing if previous scan was  
the same direction

p k WIDE

omit if same as previous scan

δx δy DXDY

" " " " "

Then

.....  
or

CENTER

After positioning PDS over appropriate  
point

.....  
or

UR (or UL, LL, LR) CORNER

After positioning PDS over appropriate  
point

.....  
or

y. x. \$CENTRE (or \$CORNER)

See [4.3]

Then loop to [8.1.2] incrementing  
n by 1 each time until all scans  
(≤ 16) have been entered

.....  
or If tape is a new tape  
Go to [8.1.3]

REWIND NTP

see [3.3] for explanation

.....  
or Otherwise position tape anywhere  
before end of double file mark



it appropriately to make the scan as if it was scanned from the UR corner. Using UL, LL or LR instead of UR would cause the final file to have the specified corner as its origin.

To transform an array scanned in y to one scanned in x

```
FIND file 1 SOURCE FIND file 2
DESTINATION TURN
```

```
*****
If file 1 is A and file 2 is B, the effect of the command TURN is,
Aij becomes Bij where subscripts represent row and column numbers.
*****
```

```
*****
Due to the simplicity of the method, the operation is prohibitively slow for
arrays >(64 by 64)
*****
```

[8.3] Tape file graphic record display.

```
*****
This utility combines SCANSALOT with the display words in SLICING (see [4.9],
[4.12]).
*****
```

454 LOAD

Loads programs

n SKIPS SET

Selects tape file, see [4.12]

FIRST

Sets scales, plots 1st line. See [4.9]

The following commands are available, and are described in [4.9]

N

Reads and plots next line

PLOT

Plots current line

AXES

BIG GRID

POINTS

FOLLOW

HISTOGRAM

n m X SIZE

n m Y SIZE

SAMPLE

Define graphic options

Plots every n<sup>th</sup> point of current line.  
n is set by n SPL : (default 5).  
like PLOT, SAMPLE does not actually  
read a record, N must be used for that  
first.

There is one additional command:

#### [ 9.4 ]

##### Coarse File

```
*****
Description: Used most often to compare scans from an electronograph with the
corresponding area on a cathode map. Disadvantage is that it is not as
accurate as ALIGN but advantage is that cathode map is only scanned once, and
file is reused for different higher resolution scans. ALIGN is generally used
to eliminate fine variations within a single scan area while the coarse-file
method is used to eliminate grosser variations across a larger area of photo-
cathode.
*****
```

```
*****
The whole or a large part of the photocathode map is scanned at a coarse
spatial resolution and stored on tape or disk. A second scan at a higher
spatial resolution is stored on disk or tape. The program creates a third file
on disk, containing photocathode map with weighted average of the 4 nearest-
neighbour points in the coarse file.
*****
```

#### [ 8.4.1 ]

##### Creating the coarse file

```
*****
Do a normal scan as in [3.2], [3.3] or [4.7.11], [4.7.2] except:
*****
ZILCH over a standard mark common to all electronographs to be scanned, (for
*****
```

```

* example, upper left fiducial mark). In addition electronographs must be
* approximately aligned with one another using the fiducial marks or photocathode
* defects. The scan must be in X and must be specified by either the CENTER or
* the UR CORNER.

```

[8.4.2] Creating the nearest-neighbour file (from the parameters of a scan, called "source scan", done from a corresponding electronograph.

```

* The source scan is done either onto tape or disk. The sample intervals must
* be less than those of the coarse file. The source scan must be in X and the
* origin must be defined by either CENTER or UR CORNER. ZIICH must be done at
* the same standard point as was done for the coarse file.

```

```

* The nearest neighbour file may now be created using the data from the coarse
* file and the parameters of the source scan.

```

```

BIGSCAN LOAD
457 LOAD
Loads coarse file routines

```

```

: If source scan on disk
FIND filename SOURCE
Selects source scan file

```

```

: Go to [8.4.3]

```

```

: Otherwise, ie. on tape
n SKIPS SET SOURCE
Selects source scan file see [4.12]
For explanation

```

```

[8.4.3] To allocate a new file in disk
LISTFILE n CREATE file 2
See [4.7.2] For explanation

```

```

FIND file 2 DESTINATION
Selects destination file

```

```

: If coarse file on disk
FIND coarse-file SECOND
Sets up coarse-file parameters

```

```


: Go to [8.4.4]

```


```

: Otherwise, ie. coarse file on
tape
n SKIPS SET SECOND
Sets up coarse file parameters

```



Example of DENSITY



Example of PLOT

[3.4.4] Then

NEAREST

Calculates and stores nearest neighbour file in file 2

.....  
\* To use file 2 created by NEAREST to take out effects of photocathode response \*  
\* from file 1 do TAYLOR LOAD PROCESSING LOAD and see section [4.8] \*  
.....

[3.5] Density histogram (Disk files only)

DHIST LOAD

Loads programs

FIND filename

Selects file

n RES

Sets bucket resolution (default 10)

DENSITY

Sets up array of 'buckets', goes through file counting number of points in file with densities in range of each bucket, and plots results as a histogram of number against density. The default bucket size n is 10 (minimum value is 1, an ordinate of 200 means that there are 200 points between the abscissa density value and that value + n

Data may be replotted as follows

a b ROI PLOT

Plots density values a to b only, without change of scale.

a b SEE PLOT

Plots density values a to b, expanding scale to fill screen.

min max Y SIZE

Sets Y scale

ENTIRE PLOT

Equivalent to 0 4095 SEE PLOT DENSITY  
Plots all values from 0 to 4095





\*\*\*\*\*  
 \* The density histogram commands documented in [ 8, 6 ] are available. \*  
 \*\*\*\*\*

.....  
 ..Either  
 .....  
 ..or  
 .....  
 ..or  
 ..a b BSET  
 .....  
 ..or  
 .....  
 ..or  
 ..a Δa LSET  
 .....  
 ..or  
 .....  
 ..or  
 ..a b XSET

Sets grey scale levels; there are 5 levels, set to be : 1)  $x \leq a$   
 2)  $a < x \leq a + \Delta a$   
 3)  $a + \Delta a \leq x \leq a + 2\Delta a$   
 4)  $a + 2\Delta a \leq x \leq b$  5)  $x \geq b$   
 where  $b = a + 3\Delta a$   
 Sets levels to be  $a + n\Delta a$ ,  $n = 0$  to  $4$   
 $\Delta a = (b - a)/5$ , and ignores points outside the range (a,b)

JPLOT  
 Clears page, plots image and puts identification on the side

GR  
 Puts grid on picture (optional)

CURSOR  
 Puts up cursor on screen:

C  
 Prints out on screen pixel Coordinate corresponding to cursor position

P  
 Types out value of array at that point

S  
 Replots at a larger scale the Segment of the picture defined by the portion of the grid it is in (you must have the grid plotted to use this).

any other key  
 Puts up the cursor again; move them a tap any key; the first and second points defined by the cursor are used as the lower, left and upper right corners of a rectangle which is then replotted on a larger scale.

[ 9. 2 ] Cursor facility

.....  
 ..Either  
 .....  
 ..or  
 .....  
 ..or  
 .....  
 ..or  
 ..any other key

DISPLAY

DISK JPLOT

Replots output on sides and grid on screen, but not picture itself. Used to select another area using the cursor without plotting the whole picture again (to save time)

This will plot the window of the file used prior to CURSOR, whereas JPLOT will just plot the latest segment or window.

[10.0]

Error conditions

If you type a command which the machine does not recognize, the machine will echo the command followed by ? and clear the push-down stack. Continue by typing in correct commands.

If you interrupt a scan to tape (by pressing Control C, then Return), the command CANCEL will move the tape back to the beginning of the file and write an extra filemark, effectively erasing the file. Note however that if you reload either BIGSCAN or TAYLOR before doing this (eg due to a system crash) the computer will not remember the current position of the tape and a CANCEL will lose the whole tape by writing a double filemark at the beginning of the tape.

The command (OFF) will send an instruction to the PDS motors cancelling any instruction to move the PDS (eg when the stage jiggles erratically after an interrupt).

1979 July 30

FILE STAR

PIXELS

256 TO 319

ROWS

320 TO 383

LEVELS

30000


24000

18000

12000

6000

0



JK



GR

FILE RICHARD

PIXELS  
228 TO 347


ROWS  
30 TO 217

LEVELS



OK

6 5/8" V SIZE PLOT



259 394 LOAD


FILE 11

PIXELS 9 TO 19


ROWS 9 TO 19

LEADS

1000  
800  
600  
400  
200  
0



OK







145 SPACE 1 DISK FLO1

COOPY  
EARTH C